

ABI-300

Manual – Single-level + Toe

simpleABI[™]

ABI-300 Unique Features:

- ✦ PC-based system
 - ✦ Easy staff training - simple to learn and use
 - ✦ Point & Click reporting right from the computer
- ✦ “Gold Standard” 8MHz Doppler for ABI indices
- ✦ “Click & Capture” PVR Waveform
- ✦ Customizable reports
- ✦ Print to any standard office printer
- ✦ Attach multiple formats directly to EMR records
- ✦ **Full 2-year warranty**



ABI-300 Kit – what’s included:

- ✦ Non-display Doppler with
 - ✦ Digital 8MHz
 - ✦ Digital PPG
- ✦ PVR waveform hardware
- ✦ Manual Aneroid
- ✦ USB drive with simpleABI software
- ✦ 6 cuffs (2x10cm, 2x12cm, 2x1.9cm)
- ✦ Mobile carry bag

ABI-300 Options:

- ✦ Computer
- ✦ DICOM Plug-in
- ✦ Additional cuff sizes

Newman  Medical

Employee Owned & Operated

<http://newman-medical.com/product-details/abi-300/>

<http://newman-medical.com/vascular-systems-resources/>

800-267-5549

WWW.NEWMAN-MEDICAL.COM