

# ABI-450CL

## Automated ABI, TBI, & Stress Testing

### ABI-450CL Highlights

- **PC-based system**
  - Easy staff training - simple to learn and use
  - Point & Click reporting right on the computer
- **Cuff-Link**
  - 8MHz Doppler
  - Automated cuff selector – no moving tubes
  - Automated inflation & deflation
- **8.5x11 Standard Reports**
  - No taping or pasting waveform strips
- **Data files easily attached to EMR records**
- **2-year warranty**

simpleABI<sup>™</sup>  
with CUFF-LINK<sup>™</sup>

### Tests Performed


- **Quick PVR \*NEW\***
  - 1-2 minute test for PAD
  - Wrap a cuff, push a button, read the results
  - No need to remove shoes or socks
- **Bi-lateral ABI**
  - Gold standard test for PAD
- **TBI**
  - Toe pressures when large arteries are incompressible
- **ABI Stress**
  - One page ABI Stress report with time graph, time to onset, maximum walking time, time to recovery

### ABI-450CL Kit – “What’s in the box”

- Cuff-Link<sup>™</sup> - 4 port cuff selector
- Cuff-Link<sup>™</sup> Remote
- simpleABI 450CL Software
- Sensitive 8MHz Doppler
- PPG Probe
- 5 cuffs (2x10cm, 2x12cm, 1x1.9cm)

### ABI-450CL Options

- Computer
- Roll Stand with Basket
- DICOM Plug-in
- Mobile Carry Bag
- Wireless Printer (wireless router needed)



<http://newman-medical.com/abi450cl.html>

<http://newman-medical.com/videos.html>

**Cuff-Link<sup>™</sup> – Push button simplicity!**