

# ABI-600CL

## Automated ABI, TBI, Segmental, & Stress Testing

### ABI-600CL Highlights

- **PC-based system**
  - Easy staff training - simple to learn and use
  - Point & Click reporting right on the computer
- **Cuff-Link**
  - 8MHz Doppler
  - Automated cuff selector – no moving tubes
  - Automated inflation & deflation
- **8.5x11 Standard Reports**
  - No taping or pasting waveform strips
- **Data files easily attached to EMR records**
- **2-year warranty**

### Tests Performed

- **Quick PVR \*NEW\***
  - 1-2 minute test for PAD
  - Wrap a cuff, push a button, read the results
  - No need to remove shoes or socks
- **Bi-lateral ABI**
  - Gold standard test for PAD
- **TBI**
  - Toe pressures when large arteries are incompressible
- **Segmental Test**
  - 3 cuff, 3 cuff with toe, & 4 cuff reports available
- **ABI Stress**
  - One page ABI Stress report with time graph, time to onset, maximum walking time, time to recovery


### ABI-600CL Kit – “What’s in the box”

- Cuff-Link™ - 8 port cuff selector
- Cuff-Link™ Remote
- simpleABI 600CL Software
- Sensitive 8MHz Doppler
- PPG Probe
- 9 cuffs (2x10cm, 4x12cm, 2x12cm long, 1x1.9cm)

### ABI-600CL Options

- Computer
- Roll Stand with Basket
- DICOM Plug-in
- Mobile Carry Bag
- Wireless Printer (wireless router needed)

simpleABI<sup>™</sup>  
with CUFF-LINK<sup>™</sup>



Single Page  
ABI Stress Report

<http://newman-medical.com/abi600cl.html>

<http://newman-medical.com/videos.html>

Newman Medical

5350 Vivian Street | Unit C | Arvada | CO | 80002

Toll-free 800-267-5549 | Fax-888-202-3627 | [www.newman-medical.com](http://www.newman-medical.com)