

# DD-700-VO

## Vascular-OB Combo Display Doppler



Newman Medical offers a full line of Doppler systems with the performance and durability you expect to meet your application - obstetrical or vascular. The **DigiDop 700 Display Digital Doppler** has unique features not found in other Doppler systems such as Digital Sound Clarity (DSC™). Each DigiDop is designed with the sensitivity you need while also designed to be tough and last through years of everyday office use. All our products are **Made in the USA in Colorado**. Each Newman Medical Doppler and probe must pass a series of comprehensive tests to assure quality and the durability of each product.



### Features:

- ✦ **3MHz & 8MHz** – Interchangeable obstetrical and vascular probes included
- ✦ **Better Sound** – Digital Sound Clarity (DSC™) produces up to 6x signal of other Dopplers
- ✦ **Better Warranty** – 5 year \$0 cost warranty – **ALL** Parts! **ALL** Labor!
- ✦ **Better Support & Service** – Knowledgeable Staff & Guaranteed Service
- ✦ Tough and durable polycarbonate and stainless steel case
- ✦ Magnetic probe holder
  - ✦ Eliminates cradle breaks
  - ✦ Reduces gel build-up and cross contamination
- ✦ User replaceable coil cord saves service downtime
- ✦ Powered by 3 AA batteries - extended battery life of 900 mins
- ✦ Dimensions (h x w x l): 150 x 65 x 35 mm (6.0 x 2.5 x 1.25 in)
- ✦ Weight: 340 grams (12 ounces)
- ✦ Shipping: 2 pounds in 12"x9"x3" box

Newman  Medical

*Employee Owned & Operated*

800-267-5549

WWW.NEWMAN-MEDICAL.COM

