


DD-901

Table Top Doppler



XD Probes provide increased sensitivity



DD-901 with Obstetrical probe



Bright LCD



Drop on recharging



**Dual Speaker
Room filling Sound**

Newman  Medical

Employee Owned & Operated

800-267-5549

WWW.NEWMAN-MEDICAL.COM



**5-year warranty covers
All parts and labor
Includes breakage!**

DD-901

Table Top Doppler



DD-901 with Vascular probes

DD-901 Features:

- ✦ Dual speaker, room-filling sound
- ✦ 5 Year warranty on all parts & labor, including breakage.
- ✦ Bright, wide angle LED backlit display
- ✦ Smart Recharge System™ (SRS)
 - ✦ Significant cost savings
 - ✦ New rechargeable batteries hold charge for years
 - ✦ May be recharged more than 2000 times
- ✦ Wall Mount and Roll Stand Options available
- ✦ Powered by 3 AA batteries
- ✦ Dimensions (h x w x l): 210 x 155 x 115 mm (8.25x 6 x 4.5 in)
- ✦ Weight with Probe: 850 grams (30 ounces)
- ✦ Shipping: 4 pounds in 13"x10"x5" box

Newman  Medical

Employee Owned & Operated

800-267-5549

WWW.NEWMAN-MEDICAL.COM

The Newman Medical **DD-901 Table Top Doppler** has breakthrough sensitivity and size. Extended Depth XD™ probes provide up to 56% increased sensitivity at shallow depths (2-5cm) where many early fetal hearts are located. Complementing the new XD probes is the state-of-the-art Digital Signal Clarity (DCS™) of the Doppler signal.

Ideal for:

- Hospitals
- Multiple practitioner offices
- Wall mounted locations



Optional Roll Stand

