

DigiDop Probes



D2

The DigiDop 2MHz probe is optimized for later term pregnancies or larger patients.



D2W

The DigiDop 2MHz waterproof probe is designed for water applications and labor deliveries.



D3

The DigiDop 3MHz probe is our general purpose obstetrical probe. It is optimized for early fetal heartbeat detection and use throughout pregnancy.



D3W

The DigiDop 3MHz waterproof probe is designed for water applications and labor deliveries.



D5

The DigiDop 5MHz probe is optimized for deeper vessels in vascular applications.



D8

The DigiDop 8MHz probe is our general purpose vascular probe. The DigiDop 8MHz probe is optimized for peripheral vessels in vascular applications.

Newman  Medical

Employee Owned & Operated

800-267-5549

WWW.NEWMAN-MEDICAL.COM



DD-PPG

The DigiDop PPG probe is for use in vascular applications. The DigiDop PPG probe is for use when a Doppler probe is insufficient – i.e. for those with incompressible arteries.

