

# DigiDop

MADE IN



THE USA



## BREAKTHROUGH

# SENSITIVITY

## XD-EXTENDED DEPTH PROBES

**XD**  
**eXtended Depth**  
**Probes**

56% more sensitive  
where you need it.

## An Honest

# 5 YEAR WARRANTY

Our 5 year warranty covers ALL labor and ALL parts Including probes and cords!




ASK YOUR SALES REPRESENTATIVE TODAY!  
FOR QUESTIONS OR SALES SUPPORT CALL

# 800-267-5549

## DD-990R Table Top

**Room Filling Sound** from the large forward-facing speaker of the DD-990R. Its bright large display is easily seen across the room.



**Large, Bright Display:** Indicates heart rate, battery status, volume setting and signal strength.

Smart  Recharge System

Our **Premium Rechargeable Batteries** stay charged when not in use and can be recharged more than 2000 times.

A **significant** cost saving and reduction of waste

## PROBE VERSATILITY: ALL PROBES CAN BE USED ON ALL OUR DOPPLERS

**OB:** Our XD probes are particularly sensitive to shallow early fetal hearts without sacrificing sensitivity at depth. Both 2 and 3 MHz are available in waterproof versions as well as standard.

**Vascular:** In addition to our sensitive and easy to hold 5 and 8 MHz probes, a unique audio PPG probe is available for detecting digital or penile flow.

MAIN UNITS	DESCRIPTION
DD-990R	DigiDop 900R Table Top Doppler
DD-300	Non-Display Digital Doppler
DD-301	Non-Display Digital Doppler with Recharger
DD-700	Display Digital Doppler
DD-701	Display Digital Doppler with Recharger

PROBES	DESCRIPTION
D3	General purpose 3MHz probe - best for early fetal hearts
D2	2MHz probe designed for deeper fetal hearts
D3W	3MHz obstetrical for water applications
D2W	2MHz obstetrical for water applications
D8	General purpose 8MHz probe - peripheral arteries
D5	5MHz probe designed for deeper vessels
DPPG	Audio PPG for TBI and other studies

## HOW TO ORDER

For product ordering code, add the probe number to the end of the main unit number with a hyphen.

For example, to order a Display Doppler (DD-700) with 8MHz probe (D8) the order code is DD-700-D8.

Visit [www.newman-medical.com](http://www.newman-medical.com) for more information and accessories such as rolling stands, carrying cases, & wall mounts