

# ABI-500CL

## Advanced – Multi-level

simpleABI<sup>™</sup>  
with CUFF:LINK<sup>™</sup>

### ABI-500CL Highlights

- ✦ **PC-based system**
  - ✦ Easy staff training - simple to learn and use
  - ✦ Point & Click reporting right on the computer
- ✦ **Cuff-Link**
  - ✦ 8MHz Doppler
  - ✦ Automated cuff selector – no moving tubes
  - ✦ Automated inflation & deflation
- ✦ **8.5x11 Standard Reports**
  - ✦ No taping or pasting waveform strips
- ✦ **Data files easily attached to EMR records**
- ✦ **2-year warranty**

### Tests Performed

- ✦ simpleABI ABI-Q<sup>™</sup> Exam
- ✦ ABI Exam
- ✦ ABI with Toe (TBI) Exam
- ✦ 3-cuff Segmental Exam
- ✦ 3-cuff with Toe Segmental Exam
- ✦ 4-cuff Segmental Exam
- ✦ 4-cuff with Toe Segmental Exam

### ABI-500CL Kit – what's included:

- ✦ simpleABI Reporting Software
- ✦ 8 port Cuff-Link<sup>™</sup> - cuff selector with remote
- ✦ Sensitive 8MHz Doppler
- ✦ PPG Probe
- ✦ 10 cuffs (2x10cm, 4x12cm, 2x12cm long, 2x1.9cm)

### ABI-500CL Options

- ✦ Computer
- ✦ DICOM Plug-in
- ✦ Mobile Carry Bag
- ✦ Additional cuff sizes



Newman  Medical

*Employee Owned & Operated*

[www.newman-medical.com/abi-500cl/](http://www.newman-medical.com/abi-500cl/)

[www.newman-medical.com/information/](http://www.newman-medical.com/information/)

800-267-5549

WWW.NEWMAN-MEDICAL.COM