

# ABI-600CL

Advanced – Multi-level + Exercise

simpleABI<sup>™</sup>  
with CUFF:LINK<sup>™</sup>

## ABI-600CL Highlights

### ✦ PC-based system

- ✦ Easy staff training - simple to learn and use
- ✦ Point & Click reporting right on the computer

### ✦ Cuff-Link

- ✦ 8MHz Doppler
- ✦ Automated cuff selector – no moving tubes
- ✦ Automated inflation & deflation

### ✦ 8.5x11 Standard Reports

- ✦ No taping or pasting waveform strips

### ✦ Data files easily attached to EMR records

### ✦ 2-year warranty

## Tests Performed

- ✦ simpleABI ABI-Q<sup>™</sup> Exam
- ✦ ABI Exam
- ✦ ABI with Toe (TBI) Exam
- ✦ 3-cuff Segmental Exam
- ✦ 3-cuff with Toe Segmental Exam
- ✦ 4-cuff Segmental Exam
- ✦ 4-cuff with Toe Segmental Exam
- ✦ ABI with Exercise (93923)
- ✦ ABI Stress-Treadmill (93924)

## ABI-600CL Kit – what's included:

- ✦ simpleABI Reporting Software
- ✦ 8-port Cuff-Link<sup>™</sup> - cuff selector with remote
- ✦ Sensitive 8MHz Doppler
- ✦ PPG Probe
- ✦ 10 cuffs (2x10cm, 4x12cm, 2x12cm long, 2x1.9cm)

## ABI-600CL Options

- ✦ Computer
- ✦ DICOM Plug-in
- ✦ Mobile Carry Bag
- ✦ Additional cuff sizes



Newman  Medical

*Employee Owned & Operated*

800-267-5549

WWW.NEWMAN-MEDICAL.COM

[www.newman-medical.com/abi-600cl/](http://www.newman-medical.com/abi-600cl/)

[www.newman-medical.com/information/](http://www.newman-medical.com/information/)