



MADE IN DENVER, CO, USA

INNOVATIVE DOPPLERS FOR TODAY'S SAVVY PHYSICIAN

- Designed, manufactured & serviced in the U.S.
- Market leading 5-year warranty
- Rechargeable systems
- Extended battery life
- Interchangeable probes



THE NEWMAN MEDICAL DIFFERENCE

QUALITY

- Unprecedented 5-year warranty
 - includes labor, parts, breakage, and probes
- All probes tested
- Acoustical output power measured for all probes
- Certified to meet international quality standard ISO 13485

SUSTAINABILITY

- Rechargeable systems include premium batteries that can be charged 2,100+ times
- Smart Recharge System charges batteries fast and won't charge if alkaline batteries are accidentally inserted

INNOVATION

- Continual product development
- Investment in doppler innovation
- Advanced features include:
 - Digital Signal Clarity
 - Wide Angle XD 3 MHz probes
 - Smart Recharge System

HISTORY

- All Dopplers are designed, manufactured, and serviced in Colorado, U.S.
- The Newman family has been involved with Doppler since it was first introduced in the 1970s



MADE IN
DENVER, CO
UNITED STATES OF AMERICA



INNOVATIVE TECHNOLOGY

Easy-to-Use | Dependable



THE DIGIDOP FAMILY OF PRODUCTS INCLUDES
COUNTLESS TOP-OF-THE-LINE FEATURES:

5-year warranty covering parts, labor, and breakage
No other Doppler manufacturer comes close!



**Large, bright display featuring
battery condition, signal
strength of Doppler or PPG
signal, and volume setting**
(some models include
calculation and display of
autocorrelated heart rate)

**Unbreakable, easy-to-clean
magnetic cradle holds probes**



**Stainless steel
microscreens protect
speakers and keep dirt out**

**Premium
batteries**

Smart  **Recharge
System**

**Digital Signal Clarity
sound processing**

**Interchangeable probes for
vascular and obstetrical use**



Professional and Sensitive Dopplers

Newman Medical



Over **45 years** in the marketplace; a testament to our longevity and success



Made in the USA for superior quality and performance



5-year all-inclusive warranty for maximum protection and minimized downtime, headaches



Best-in-class customer support; we always answer the phone

DigiDop Classic

- Standard size
- Superior battery endurance compared to competitors



DigiDop II

- Compact size
- Comfortable palm-fitting design
- Equivalent battery endurance to competitors



DigiDop 901

- Ideal for hospital use
- Clear stereo sound from dual top speakers
- Two probe cradles
- Standard drop-on charging pyramid
- Optional wall mount recharger



Vascular Probes



Model D8
8MHz Vascular Probe
or
Model D5
5MHz Vascular Probe



Model D8-SD
8 MHz SteadyDop Vascular Probe
or
Model D5-SD
5 MHz SteadyDop Vascular Probe



Model DPPG
Audio TBI PPG Probe

Fetal Probes



Model D2
2MHz Fetal Probe
or
Model D3
3MHz Fetal Probe



Model D2W
2MHz Waterproof Fetal Probe
or
Model D3W
3MHz Waterproof Fetal Probe

Interchangeable probes can be used on any DigiDop

Vascular



- Great sensitivity
- Ergonomic shape easy to hold, place in sterile glove finger for use in sterile fields
 - **D8** 8 MHz - best for shallow vessels like radial, ankle, foot, digits, penile
 - **D8-SD** - 8 MHz SteadyDop - best for shallow vessels like radial, ankle, foot, digits, penile
 - **D5** - 5 MHz - wider beam, deeper penetration. Good for brachial, carotid, femoral, deep veins
 - **D5-SD** - 5 MHz SteadyDop - wider beam, deeper penetration. Good for brachial, carotid, femoral, deep veins

Fetal



- Great sensitivity
- Smooth shape, comfortable for patient, easy to hold and clean
 - **D3** - 3 MHz - wide angle, ideal for early fetal hearts
 - **D3W** - a unique 3 MHz waterproof probe
 - **D2** - 2MHz - deeper penetration
 - **D2W** - 2 MHz - waterproof

DPPG



- Unique audio PPG probe for use on digits. Great for toe-brachial studies

How to build your Doppler

DOPPLER CODE

PROBE CODE

Doppler

DD-300	Basic
DD-301	Basic + Recharge
DD-700	HR Display
DD-701	HR Display + Recharge
DD-901	Tabletop + HR Display + Recharge

DigiDop II

DD-330	Basic
DD-330R	Basic + Recharge
DD-770	HR Display
DD-770R	HR Display + Recharge

Probe

D2	2MHz fetal probe
D2W	2MHz waterproof fetal probe
D3	3MHz fetal probe
D3W	3MHz waterproof fetal probe
D5	5MHz vascular probe
D5-SD	5MHz SteadyDop vascular probe
D8	8MHz vascular probe
D8-SD	8MHz SteadyDop vascular probe
DPPG	Audio PPG probe - includes toe cuff

DigiDop 2-in-1 Dopplers

with interchangeable fetal & vascular probes

Newman Medical

DigiDop

DigiDop Tabletop



(Larger size | Longer battery life)

DigiDop "Classic"



(Standard size | Longer battery life)

DigiDop II



(Pocket size | Standard battery life)

DigiDop Tabletop

DigiDop "Classic"

DigiDop II

DD-901

DD-300

DD-301

DD-700

DD-701

DD-330

DD-330R

DD-770

DD-770R

Heartrate Display



Non-Display



Rechargeable



Alkaline



APPLICATION

SKU

DESCRIPTION



D2

2 MHz fetal probe

D2W

2 MHz waterproof fetal probe

D3

3 MHz fetal probe

D3W

3 MHz waterproof fetal probe



D5

5 MHz vascular probe

D5-SD

5 MHz SteadyDop vascular probe

D8

5 MHz vascular probe

D8-SD

8 MHz SteadyDop vascular probe

DPPG

Audio PPG probe - includes toe cuff

NEW

INTRODUCING THE STEADYDOP VASCULAR PROBE

Easy ABI Measurements, Even for Beginners:

The SteadyDop probe eliminates movement errors and angle guesswork, making accurate vascular assessments simple with your simpleABI system* – no prior Doppler experience necessary.

*The SteadyDop vascular probe is available as a separate accessory for use with simpleABI PAD testing systems and DigiDop Doppler units.



POPULAR ACCESSORIES

For a complete listing of accessories go to
newman-medical.com or call customer service

MADE IN

USA



WALL MOUNT



ROLL STANDS



ABI-TBI KITS

Newman  Medical

MADE IN
DENVER, CO
UNITED STATES OF AMERICA

DROP-IN
HOLDER/
CHARGER/

Smart  Recharge
System



NM0107