




Primary Care



Podiatry



Wound Care



Community Health

CLINIC PLACEMENT



Rapid ABI-Q

EXAMS PERFORMED



American-Made Quality



Includes labor, parts, and accidental breakage



ABI-Q Rapid Screen



DICOM Compatible

ABI-Q

The simpleABI-Q is the fastest path from patient encounter to PAD answer – no Doppler required. Designed specifically for the patients most likely to be missed, it uses PVR waveform technology to deliver accurate results even in diabetics and patients with calcified arteries where traditional Doppler ABI falls short. Portable, staff-friendly, and reimbursable, the ABI-Q fits into any clinical environment – from primary care to mobile screenings – and adds exercise testing capability without adding complexity. When speed and accessibility matter most, this is the system that removes every barrier between your at-risk patients and an answer.

ABI-Q Includes:

- simpleABI software
- PVR
- Cuffs
 - 2 x 10cm EZ-Clean Vasocuffs
- Inflation bulb
- Carry bag
- User manual / Quick reference guide
- Two-Year Warranty
 - Includes all parts, labor, and accidental breakage

ABI-Q Exams

- ABI-Q Rapid Screen



simpleABI Vascular Systems: Built for the Patients Who Can't Wait

Proudly Made in the USA

- Precision-engineered in Colorado, not outsourced, not compromised.
- Rugged enough for the demands of real clinical environments.
- 2-year full-coverage warranty, because reliability shouldn't be optional.

Newman Medical: 50+ Years of Vascular Expertise

- Pioneering ultrasound diagnostics since the 1970s.
- Innovation driven by one goal: better outcomes for at-risk patients.
- Long-term partnerships built on results, not promises.

A Support Team as Committed as You Are

- Deep clinical expertise on call when questions arise.
- Field representatives dedicated to your practice's success.
- Robust online and in-person resources, always within reach.



ABI-Q



Exams Performed



ABI-Q Rapid Screen



Single-Level PVR

Ordering Information

System

| SKU | DESCRIPTION |
|-------|---------------------------------|
| ABI-Q | simpleABI-Q Rapid Screen System |

Accessories

| SKU | DESCRIPTION |
|----------|---|
| ABI-120 | Microsoft® Windows® PC for use with simpleABI Systems |
| ACC-290 | Carry Bag for simpleABI-Q Systems |
| STND-130 | simpleABI Roll Stand for simpleABI-Q |
| ACC-230 | Memory Stick, USB with simpleABI Software |
| PUMP-100 | Manual cuff inflator |



Replacement Cuffs

| SKU | DESCRIPTION |
|--------------|--|
| CUFF-200 | EZ Clean Vascular Cuff, 10cm (ankle, small arm) |
| CUFF-210 | EZ Clean Vascular Cuff, 12cm (arm, ankle, calf) |
| CUFF-SET-242 | EZ Clean Vascular Cuff Set - 2x10cm, 4x12cm, 2x12cm - 8/pk |

ABI-Q Technical Specifications

| | |
|--|---|
| Weight: | 3.6 kgs (8 pounds) |
| Level of Protection against electrical shock: | Type B Applied Part, Class II Equipment |
| Comply with the following Standards: | IEC60101-1, IEC60601-2 |
| Transport/Storage Temperature: | -20° ~ 50°C (-4° ~ 122°F) |
| Transport/Storage Humidity: | 5% ~ 90%, non-condensing |
| GTIN: | 00852785008661 |

REIMBURSEMENT INFORMATION

Newman Medical's simpleABI vascular systems frequently utilize CPT codes 93922, 93923, and 93924 for billing and reimbursement purposes.

To obtain expert guidance on the appropriateness of these codes for specific systems, we highly recommend consulting with one of our experienced vascular specialists.

