



Primary Care



Podiatry



Wound Care



Community Health

CLINIC PLACEMENT



Rapid ABI-Q

EXAMS PERFORMED



American-Made Quality



Includes labor, parts, and accidental breakage



ABI-Q Rapid Screen



DICOM Compatible

ABI-Q

The ABI-Q Rapid-Test System delivers fast, accurate PAD diagnosis in seconds, right in your office. Get comprehensive reports quickly, enabling early detection and improved patient outcomes. Cost-effective and portable, ABI-Q is ideal for primary care settings, especially for patients with diabetes or those with calcified arteries. Simplify testing, enhance care, and benefit from reimbursement opportunities with ABI-Q.

ABI-Q Exams

- ABI-Q Rapid Screen

ABI-Q Includes:

- simpleABI software
- PVR
- Cuffs
 - 2 x 10cm EZ-Clean Vasocuffs
- Inflation bulb
- Carry bag
- User manual / Quick reference guide
- Two-Year Warranty
 - Includes all parts, labor, and accidental breakage



simpleABI Vascular Systems: Quality You Can Trust

Made in the USA

- High-quality components, manufactured in Colorado.
- Durable design for busy clinical settings.
- 2-year warranty with complete coverage.




Newman Medical: Your Partner in Vascular Care

- Leaders in ultrasound diagnostics since the 1970s.
- Committed to innovation and providing exceptional vascular products.
- Focused on service and long-term customer relationships.

First-Class Support

- Experienced customer service team with extensive industry knowledge.
- Dedicated field representatives for personalized assistance.
- Comprehensive clinical support resources, both online and in-person.


ABI-Q



Exams Performed

5 min

ABI-Q Rapid Screen



Single-Level PVR

Ordering Information

System

SKU

ABI-Q

DESCRIPTION

simpleABI-Q Rapid Screen System

Accessories

SKU

ABI-120

ACC-290

STND-130

ACC-230

PUMP-100

DESCRIPTION

Microsoft® Windows® PC for use with simpleABI Systems

Carry Bag for simpleABI-Q Systems

simpleABI Roll Stand for simpleABI-Q

Memory Stick, USB with simpleABI Software

Manual cuff inflator

Replacement Cuffs

SKU

CUFF-200

CUFF-210

CUFF-SET-242

DESCRIPTION

EZ Clean Vascular Cuff, 10cm (ankle, small arm)

EZ Clean Vascular Cuff, 12cm (arm, ankle, calf)

EZ Clean Vascular Cuff Set - 2x10cm, 4x12cm, 2x12cm - 8/pk



ABI-Q Technical Specifications

Weight:

3.6 kgs (8 pounds)

Level of Protection against electrical shock:

Type B Applied Part, Class II Equipment

Comply with the following Standards:

IEC60101-1, IEC60601-2

Transport/Storage Temperature:

-20° ~ 50°C (-4° ~ 122°F)

Transport/Storage Humidity:

5% ~ 90%, non-condensing

GTIN:

00852785008661

REIMBURSEMENT INFORMATION

Newman Medical's simpleABI vascular systems frequently utilize CPT codes 93922, 93923, and 93924 for billing and reimbursement purposes.

To obtain expert guidance on the appropriateness of these codes for specific systems, we highly recommend consulting with one of our experienced vascular specialists.

