



DICOM for simpleABI Overview

DICOM for simpleABI is an application that helps you send documents or images directly to PACS. DICOM for simpleABI can combine a multitude of actions into powerful DICOM workflows.



DICOM for simpleABI Unique Features

- **Programmable background operation**
- **Parses documents and queries worklist and stores to PACS**
- **Hold and Re attempt**
 - DICOM for simpleABI can wait and intermittently requery until associated exam records and images arrive in PACS, so it's perfect for attaching referral documents to future diagnostic test results.
- **Non intrusive and reliable**
 - DICOM for simpleABI embodies a unique approach to workflows by providing robust "fire and forget" DICOM functionality to applications and devices - no more waiting until things get stored.
- **Standard DICOM formats or Structured Reporting**
 - DICOM for simpleABI is a configurable solution that can export exam information in standard DICOM or configurable structured report format.
- **Site Licenses Available**
 - For heavy users of DICOM with multiple units, site licenses are available.
 - Up to 25 installs for the price of 10.
 - Up to 100 installs for the price of 25.



Order Information:

SKU: SFTW-100

Description: DICOM Software Integration.
Full worklist capability. Image & Structured Reporting

