# simpleABI-400CL Automated Cuff-Link System

















**Primary Care** 

**Podiatry** 

Wound Care Community Health

CLINIC PLACEMENT

**EXAMS PERFORMED** 





Includes labor, parts, and accidental breakage



**Automatic Cuff Selector** 



**Push-Button Remote** 



simpleABI Roll Stand



ABI-Q Rapid Screen



The simpleABI-400CL Automated Cuff-Link System offers accurate and efficient PAD testing for confident, early diagnosis. With versatile testing options, streamlined workflows, and simplified reimbursement processes, simpleABI enhances convenience and safety for both clinicians and patients. Backed by a 2-year warranty, it's the smart choice for proactive PAD management and improved patient outcomes.

# ABI-400CL Exams

- ABI-Q Rapid Screen
- Single-Level
- TBI with PPG

# ABI-400CL Includes:

- simpleABI software
- 4-port Cuff-Link with remote
- 8 MHz Doppler probe
- DPPG probe
- Cuffs
  - o 2 x 10cm EZ-Clean Vasocuffs
  - 2 x 12cm EZ-Clean Vasocuffs 2 x 1.9cm disposable digit cuffs
- Adjustable-height roll stand
- User manual / Quick reference guide
- Two-Year Warranty
  - Includes all parts, labor, and accidental breakage



# simpleABI Vascular Systems: Quality You Can Trust

# Made in the USA

- High-quality components, manufactured in Colorado.
- Durable design for busy clinical settings.
- 2-year warranty with complete coverage.

# **Newman Medical: Your Partner in Vascular Care**

- Leaders in ultrasound diagnostics since the 1970s.
- Committed to innovation and providing exceptional vascular products.
- Focused on service and long-term customer relationships.

# **First-Class Support**

- Experienced customer service team with extensive industry knowledge.
- Dedicated field representatives for personalized assistance.
- Comprehensive clinical support resources, both online and in-person.

# SIMPLIFY VASCULAR ASSESSMENTS WITH THE STEADYDOP PROBE



Its unique stethoscope-like grip and optimized design ensure easy vessel tracking. Compatible with Newman Medical's DigiDop & simpleABI systems. Optional accessory. Contact your Newman Medical rep for details.


# simpleABI-400CL Automated Cuff-Link System Newman ) Medical



# ABI-400CL



**Exams Performed** 





# **Ordering Information**

System

SKU **DESCRIPTION** 

ABI-400CL simpleABI Single-Level Cuff-Link System

**Cuff-Link Accessory Options** 

**DESCRIPTION** SKU

ABI-120 Microsoft® Windows® PC for use with simpleABI Systems

STND-500 Variable-Height Roll Stand for Cuff-Link Systems

Large simpleABI carry bag ACC-190 D5 5 Mhz vascular probe only

D5-SD 5 MHz SteadyDop vascular probe only

D8 8 MHz vascular probe only

D8-SD 8 MHz SteadyDop vascular probe only

SPPG PPG Probe for use with Cuff-Link Systems - includes toe cuff

SFTW-100 DICOM software integration

ARI-WARR 1 year additional warranty for simpleABI systems

# **Replacement Cuffs**

**DESCRIPTION** 

CUFF-200 EZ Clean Vascular Cuff, 10cm (ankle, small arm) CUFF-210 EZ Clean Vascular Cuff, 12cm (arm, ankle, calf) CUFF-210L EZ Clean Vascular Cuff, 12cm long (above knee)

CUFF-220D Digit Cuff, 2.5cm - Latex-Free CUFF-240D Digit Cuff, 1.9cm - Latex-Free

CUFF-SET-242 EZ Clean Vascular Cuff Set - 2x10cm, 4x12cm, 2x12cm - 8/pk



# SIMPLIFY VASCULAR ASSESSMENTS WITH THE STEADYDOP PROBE



Its unique stethoscope-like grip and optimized design ensure easy vessel tracking. Compatible with Newman Medical's DigiDop & simpleABI systems. Optional accessory. Contact your Newman Medical rep for details.

# simpleABI Cuff-Link Technical Specifications

Height: 107 cm (42 inches) Weight: 5.5 kgs (12 pounds)

Level of Protection against electrical shock: Type B Applied Part, Class II Equipment Comply with the following Standards: IEC60101-1, IEC60601-2, IEC60601-2-37

Operating Temperature:  $10^{\circ} \sim 40^{\circ}\text{C} (50^{\circ} \sim 104^{\circ}\text{F})$ 

**Operating Humidity:**  $30\% \sim 75\%$ 

**Transport/Storage Temperature:** -20° ~ 50°C (-4° ~ 122°F) **Transport/Storage Humidity:** 5% ~ 90%, non-condensing

Battery Life: 10 hours

### REIMBURSEMENT INFORMATION

Newman Medical's simpleABI vascular systems frequently utilize CPT codes 93922, 93923, and 93924 for billing and reimbursement purposes.

To obtain expert guidance on the appropriateness of these codes for specific systems, we highly recommend consulting with one of our experienced vascular specialists.

