

Primary Care Podiatry Wound Care Community Health

CLINIC PLACEMENT



Rapid ABI-Q Single-Level TBI w/PPG

EXAMS PERFORMED



American-Made Quality



Includes labor, parts, and accidental breakage



Automatic Cuff Selector



Push-Button Remote



simpleABI Roll Stand



ABI-Q Rapid Screen



DICOM Compatible

## simpleABI with CUFF-LINK

The simpleABI-400CL is the entry point into automated Cuff-Link PAD testing – and it's a substantial one. Combining the ABI-Q rapid Doppler-free screen with single-level Doppler ABI, TBI with waveform output, and exercise testing, it covers the core clinical scenarios that drive the majority of PAD diagnoses in primary and specialty care.

The automated cuff sequencing and push-button remote mean any staff member can run a complete exam with confidence, and the included PC and roll stand mean the system is ready to work the moment it arrives. For practices building a vascular testing capability from the ground up, the simpleABI-400CL is where comprehensive PAD care begins.

### ABI-400CL Includes:

- simpleABI software
- 4-port Cuff-Link with remote
- 8 MHz Doppler probe
- DPPG probe
- Cuffs
  - 2 x 10cm EZ-Clean Vasocuffs
  - 2 x 12cm EZ-Clean Vasocuffs
  - 2 x 1.9cm disposable digit cuffs
- Adjustable-height roll stand
- User manual / Quick reference guide
- Two-Year Warranty
  - Includes all parts, labor, and accidental breakage

### ABI-400CL Exams

- ABI-Q Rapid Screen
- Single-Level
- TBI with PPG



## simpleABI Vascular Systems: Built for the Patients Who Can't Wait

### Proudly Made in the USA

- Precision-engineered in Colorado, not outsourced, not compromised.
- Rugged enough for the demands of real clinical environments.
- 2-year full-coverage warranty, because reliability shouldn't be optional.

### Newman Medical: 50+ Years of Vascular Expertise

- Pioneering ultrasound diagnostics since the 1970s.
- Innovation driven by one goal: better outcomes for at-risk patients.
- Long-term partnerships built on results, not promises.

### A Support Team as Committed as You Are

- Deep clinical expertise on call when questions arise.
- Field representatives dedicated to your practice's success.
- Robust online and in-person resources, always within reach.




## ABI-400CL



### Exams Performed



ABI-Q Rapid Screen



Single-Level



TBI with PPG

## Ordering Information

### System

SKU	DESCRIPTION
ABI-400CL	simpleABI Single-Level Cuff-Link System

### Cuff-Link Accessory Options

SKU	DESCRIPTION
ABI-120	Microsoft® Windows® PC for use with simpleABI Systems
STND-500	Variable-Height Roll Stand for Cuff-Link Systems
ACC-190	Large simpleABI carry bag
D5	5 Mhz vascular probe only
D5-SD	5 MHz SteadyDop vascular probe only
D8	8 MHz vascular probe only
D8-SD	8 MHz SteadyDop vascular probe only
SPPG	PPG Probe for use with Cuff-Link Systems - includes toe cuff
SFTW-100	DICOM software integration
ABI-WARR	1 year additional warranty for simpleABI systems

### Replacement Cuffs

SKU	DESCRIPTION
CUFF-200	EZ Clean Vascular Cuff, 10cm (ankle, small arm)
CUFF-210	EZ Clean Vascular Cuff, 12cm (arm, ankle, calf)
CUFF-210L	EZ Clean Vascular Cuff, 12cm long (above knee)
CUFF-220D	Digit Cuff, 2.5cm - Latex-Free
CUFF-240D	Digit Cuff, 1.9cm - Latex-Free
CUFF-SET-242	EZ Clean Vascular Cuff Set - 2x10cm, 4x12cm, 2x12cm - 8/pk



## SIMPLIFY VASCULAR ASSESSMENTS WITH THE STEADYDOP PROBE



Its unique stethoscope-like grip and optimized design ensure easy vessel tracking. Compatible with Newman Medical's DigiDop & simpleABI systems. Optional accessory. Contact your Newman Medical rep for details.



## simpleABI Cuff-Link Technical Specifications

<b>Height:</b>	107 cm (42 inches)
<b>Weight:</b>	5.5 kgs (12 pounds)
<b>Level of Protection against electrical shock:</b>	Type B Applied Part, Class II Equipment
<b>Comply with the following Standards:</b>	IEC60101-1, IEC60601-2, IEC60601-2-37
<b>Operating Temperature:</b>	10° ~ 40°C (50° ~ 104°F)
<b>Operating Humidity:</b>	30% ~ 75%
<b>Transport/Storage Temperature:</b>	-20° ~ 50°C (-4° ~ 122°F)
<b>Transport/Storage Humidity:</b>	5% ~ 90%, non-condensing
<b>Battery Life:</b>	10 hours

## REIMBURSEMENT INFORMATION

Newman Medical's simpleABI vascular systems frequently utilize CPT codes 93922, 93923, and 93924 for billing and reimbursement purposes.

To obtain expert guidance on the appropriateness of these codes for specific systems, we highly recommend consulting with one of our experienced vascular specialists.

