
simpleABI-600CL Automated Cuff-Link System



















PLACEMENT

EXAMS PERFORMED



American-Made Quality



Includes labor, parts, and accidental breakage



Automatic Cuff Selector



Push-Button Remote



simpleABI Roll Stand



ABI-Q Rapid Screen



The simpleABI-600CL Automated Cuff-Link System offers accurate and efficient PAD testing for confident, early diagnosis. With versatile testing options, streamlined workflows, and simplified reimbursement processes, simpleABI enhances convenience and safety for both clinicians and patients. Backed by a 2-year warranty, it's the smart choice for proactive PAD management and improved patient outcomes.

ABI-600CL Exams

- ABI-Q Rapid Screen
- Single-Level
- Multi-Level
- TBI with PPG
- Exercise & Stress

ABI-600CL Includes:

- simpleABI software
- 4-port Cuff-Link with remote
- 8 MHz Doppler probe
- DPPG probe
- Cuffs

 - 2 x 10cm EZ-Clean Vasocuffs 4 x 12cm EZ-Clean Vasocuffs 2 x 12cm long EZ-Clean Vasocuffs
 - 2 x 1.9cm disposable digit cuffs
- Adjustable-height roll stand User manual / Quick reference guide
- Two-Year Warranty
 - Includes all parts, labor, and accidental breakage

simpleABI Vascular Systems: Quality You Can Trust

Made in the USA

- High-quality components, manufactured in Colorado.
- Durable design for busy clinical settings.
- 2-year warranty with complete coverage.

Newman Medical: Your Partner in Vascular Care

- Leaders in ultrasound diagnostics since the 1970s.
- Committed to innovation and providing exceptional vascular products.
- Focused on service and long-term customer relationships.

First-Class Support

- Experienced customer service team with extensive industry knowledge.
- Dedicated field representatives for personalized assistance.
- Comprehensive clinical support resources, both online and in-person.

SIMPLIFY VASCULAR ASSESSMENTS WITH THE STEADYDOP PROBE



Its unique stethoscope-like grip and optimized design ensure easy vessel tracking. Compatible with Newman Medical's DigiDop & simpleABI systems. Optional accessory. Contact your Newman Medical rep for details.





simpleABI-600CL Automated Cuff-Link System Newman) Medical















Ordering Information

System

SKU **DESCRIPTION**

ABI-600CL simpleABI Multi-Level with Exercise Cuff-Link System

Cuff-Link Accessory Options

SKU **DESCRIPTION**

Microsoft® Windows® PC for use with simpleABI Systems ABI-120

STND-500 Variable-Height Roll Stand for Cuff-Link Systems

ACC-190 Large simpleABI carry bag D5 5 Mhz vascular probe only

D5-SD 5 MHz SteadyDop vascular probe only

D8 8 MHz vascular probe only

D8-SD 8 MHz SteadyDop vascular probe only

SPPG PPG Probe for use with Cuff-Link Systems - includes toe cuff

SFTW-100 DICOM software integration

ABI-WARR 1 year additional warranty for simpleABI systems

Replacement Cuffs

SKU **DESCRIPTION** EZ Clean Vascular Cuff, 10cm (ankle, small arm) CUFF-200 CUFF-210 EZ Clean Vascular Cuff, 12cm (arm, ankle, calf) CUFF-210L EZ Clean Vascular Cuff, 12cm long (above knee) EZ-Clean Vascular Cuff, 17cm Contoured (thigh) CUFF-230 CUFF-250 EZ-Clean Vascular Cuff, 7cm long (metatarsal)

Digit Cuff, 2.5cm - Latex-Free CUFF-220D CUFF-240D Digit Cuff, 1.9cm - Latex-Free

CUFF-SET-242 EZ Clean Vascular Cuff Set - 2x10cm, 4x12cm, 2x12cm - 8/pk



SIMPLIFY VASCULAR ASSESSMENTS WITH THE STEADYDOP PROBE



Its unique stethoscope-like grip and optimized design ensure easy vessel tracking. Compatible with Newman Medical's DigiDop & simpleABI systems. Optional accessory. Contact your Newman Medical rep for details.

ABI-600CL Technical Specifications

Height: 107 cm (42 inches) Weight: 5.5 kgs (12 pounds)

Level of Protection against electrical shock: Type B Applied Part, Class II Equipment

Comply with the following Standards: IEC60101-1, IEC60601-2, IEC60601-2-37

Operating Temperature: $10^{\circ} \sim 40^{\circ} \text{C} (50^{\circ} \sim 104^{\circ} \text{F})$

Operating Humidity: $30\% \sim 75\%$

Transport/Storage Temperature: $-20^{\circ} \sim 50^{\circ}\text{C} (-4^{\circ} \sim 122^{\circ}\text{F})$ Transport/Storage Humidity: 5% ~ 90%, non-condensing

> Battery Life: 10 hours GTIN: 852785008609

REIMBURSEMENT INFORMATION

Newman Medical's simpleABI vascular systems frequently utilize CPT codes 93922, 93923, and 93924 for billing and reimbursement purposes.

To obtain expert guidance on the appropriateness of these codes for specific systems, we highly recommend consulting with one of our experienced vascular specialists.

