simpleABI-300 Handheld PVR System

















PLACEMENT CLINIC

EXAMS PERFORMED



American-Made Quality





simple ABI

The simpleABI Handheld PVR system is a portable, clinically precise, and cost-effective solution for detecting Peripheral Artery Disease (P.A.D.). It offers enhanced manual control and clinical accuracy, simplifies reimbursement, and comes with a 2-year warranty.

The small design of simpleABI-300 makes it suitable for various healthcare practices, including mobile imaging services and community health fairs. The system is proudly made in the USA and meets the highest-quality standards.

ABI-300 Exams

- Single-Level
- TBI with PPG

ABI-300 Includes:

- Cuffs
 - 2 x 10cm EZ-Clean Vasocuffs
 - o 2 x 12cm EZ-Clean Vasocuffs
 - o 2 x 2.5cm disposable digit cuffs
- Two-Year Warranty
 - o Includes all parts, labor, and accidental breakage

American-Made Quality

- simpleABI systems are manufactured utilizing the finest quality components and are produced outside of Denver, Colorado.
- Our products are designed to endure years of use, even in the busiest of clinical settings.
- Our industry-leading 2-year warranty provides comprehensive coverage, including labor, parts, and accidental breakage, ensuring complete peace of mind.

Newman Medical History & Commitment

- Since the 1970s, the Newman family has been at the forefront of ultrasound diagnostics.
- Newman Medical is committed to ongoing product development, guaranteeing the consistent availability of top-notch vascular diagnostic products in the market.
- Newman Medical prioritizes customer service and endeavors to establish lasting partnerships with its clientele.

First-Class Support

- Newman Medical takes great pride in its customer service team, which comprises some of the most seasoned professionals in the industry. On average, our team has over 25 years of experience.
- Our highly proficient field representatives are always ready to offer physicians exceptional support.
- Newman Medical provides physicians with a range of clinical support options, including online and in-person resources, to ensure comprehensive assistance.

simpleABI-300 Handheld PVR System



ABI-300



Exams Performed





Ordering Information

System

ABI-300 simpleABI Single-Level Cuff-Link System with TBI

Accessories

DESCRIPTION SKU

ABI-120 Microsoft® Windows® PC for use with simpleABI Systems

STND-130 simpleABI Roll Stand for Manual Systems

simpleABI-300 carry bag ACC-175 PUMP-100 Manual cuff inflator

Replacement Cuffs

DESCRIPTION

CUFF-200 EZ Clean Vascular Cuff, 10cm (ankle, small arm) CUFF-210 EZ Clean Vascular Cuff, 12cm (arm, ankle, calf) CUFF-210L EZ Clean Vascular Cuff, 12cm long (above knee)

Digit Cuff, 2.5cm - Latex-Free CUFF-220D Digit Cuff, 1.9cm - Latex-Free CUFF-240D

CUFF-SET-242 EZ Clean Vascular Cuff Set - 2x10cm, 4x12cm, 2x12cm - 8/pk



Weight: 3.6 kgs (8 pounds)

Level of Protection against electrical shock: Type B Applied Part, Class II Equipment Comply with the following Standards: IEC60101-1, IEC60601-2, IEC60601-2-37

Operating Temperature: $10^{\circ} \sim 40^{\circ}\text{C} (50^{\circ} \sim 104^{\circ}\text{F})$

Operating Humidity: 30% ~ 75%

Transport/Storage Temperature: -20° ~ 50°C (-4° ~ 122°F) Transport/Storage Humidity: 5% ~ 90%, non-condensing

Battery Life: 800, 1 minute exams GTIN: 00852785008647

REIMBURSEMENT INFORMATION

Newman Medical's simpleABI vascular systems frequently utilize CPT codes 93922, 93923, and 93924 for billing and reimbursement purposes.

To obtain expert guidance on the appropriateness of these codes for specific systems, we highly recommend consulting with one of our experienced vascular specialists.

