simpleABI-600CL Automated Cuff-Link System Newman Medical Specialists Cardiology Vascular Surgery Interventional Radiology Vascular Surgery Interventional Radiology Rapid ABI-Q Single-Level Multi-Level TBI w/PPG Exercise

EXAMS

PERFORMED

simple AB

CLINIC PLACEMENT



American-Made Quality



Includes labor, parts, and accidental breakage



Automatic Cuff Selector



Push-Button Remote





ABI-Q Rapid Screen



simple ABI

The simpleABI Automated Cuff-Link System provides tailored testing capabilities and clinical precision for early detection of Peripheral Artery Disease (P.A.D.). It offers versatile testing options, efficiency, and simplifies reimbursement processes, enhancing convenience and safety. With a 2-year warranty, it is the ultimate choice for proactive P.A.D. diagnosis and patient well-being.

ABI-600CL Exams

- ABI-Q Rapid Screen
- Single-Level
- Multi-Level
- TBI with PPG
- Exercise & Stress

ABI-600CL Includes:

- Cuffs
 - 2 x 10cm EZ-Clean Vasocuffs
 - 4 x 12cm EZ-Clean Vasocuffs
 - 2 x 12cm EZ
 - Clean Vasocuffs
 2 x 2 Fore dispersible digit suffs
 - 2 x 2.5cm disposable digit cuffs
- Two-Year Warranty
- Includes all parts, labor, and accidental breakage

American-Made Quality

- simpleABI systems are manufactured utilizing the finest quality components and are produced outside of Denver, Colorado.
- Our products are designed to endure years of use, even in the busiest of clinical settings.
- Our industry-leading 2-year warranty provides comprehensive coverage, including labor, parts, and accidental breakage, ensuring complete peace of mind.

Newman Medical History & Commitment

- Since the 1970s, the Newman family has been at the forefront of ultrasound diagnostics.
- Newman Medical is committed to ongoing product development, guaranteeing the consistent availability of high-quality vascular diagnostic products in the market.
- Newman Medical has a relentless focus on service and establishing long term relationships with our customers

First-Class Support

- Newman Medical takes great pride in its customer service team, which comprises some of the most seasoned professionals in the industry. On average, our team has over 25 years of experience.
- Our highly proficient field representatives are always ready to offer physicians exceptional support.
 Newman Medical provides physicians with a range of clinical support options, including online and in-person resources, to ensure comprehensive assistance.
- INTOOD4-F

simpleABI-600CL Automated Cuff-Link System Newman) Medical

CUFF-250

CUFF-220D


CUFF-240D

CUFF-SET-242



ABI-Q Rapid Screen









	Ordering Information	
System <u>SKU</u> ABI-600CL	DESCRIPTION simpleABI Multi-Level with Exercise Cuff-Link System	
Accessories <u>SKU</u> ABI-120 STND-500 ACC-190 D5 SPPG SFTW-100	DESCRIPTION Microsoft® Windows® PC for use with simpleABI Systems Variable-Height Roll Stand for Cuff-Link Systems Large simpleABI carry bag 5Mhz vascular probe only PPG Probe for use with Cuff-Link Systems - includes toe cuff DICOM software integration	
ABI-WARR Replacement Cuffs SKU CUFF-200 CUFF-210 CUFF-210L CUFF-230	1 year additional warranty for simpleABI systems DESCRIPTION EZ Clean Vascular Cuff, 10cm (ankle, small arm) EZ Clean Vascular Cuff, 12cm (arm, ankle, calf) EZ Clean Vascular Cuff, 12cm long (above knee) EZ-Clean Vascular Cuff, 17cm Contoured (thigh)	

EZ-Clean Vascular Cuff, 7cm long (metatarsal)

EZ Clean Vascular Cuff Set - 2x10cm, 4x12cm, 2x12cm - 8/pk

Digit Cuff, 2.5cm - Latex-Free

Digit Cuff, 1.9cm - Latex-Free

ABI-600CL Technical Specifications

Level of Protection against electrical shock: Type B Applied Part, Class II Equipment **Operating Humidity:** 30% ~ 75% Transport/Storage Temperature: -20° ~ 50°C (-4° ~ 122°F) Transport/Storage Humidity: 5% ~ 90%, non-condensing Battery Life: 10 hours

Height: 107 cm (42 inches) Weight: 5.5 kgs (12 pounds) Comply with the following Standards: IEC60101-1, IEC60601-2, IEC60601-2-37 **Operating Temperature:** $10^{\circ} \sim 40^{\circ}C (50^{\circ} \sim 104^{\circ}F)$ GTIN: 852785008609

REIMBURSEMENT INFORMATION

Newman Medical's simpleABI vascular systems frequently utilize CPT codes 93922, 93923, and 93924 for billing and reimbursement purposes.

To obtain expert guidance on the appropriateness of these codes for specific systems, we highly recommend consulting with one of our experienced vascular specialists.

