

Discover the DigiDop Doppler that closely matches your Huntleigh Doppler below.



Includes labor, parts, and accidental breakage



American-Made Quality



Includes cords and probes



Waterproof Option



Rechargeable Option

VASCULAR DOPPLERS

HUNTLEIGH Model	REPLACEMENT Newman Medical DigiDop Classic	dopplex	REPLACEMENT Newman Medical DigiDop II
NON DISPLAY D900 Dopplex D900	DD-300	or	DD-330
SD2 Dopplex SD2	DD-700	or	DD-770
DISPLAY DMX	DD-700	or	DD-770
DMXR Dopplex DMX/DMXR	DD-701	or	DD-770R

FETAL DOPPLERS

HUNTLEIGH Model	REPLACEMENT Newman Medical DigiDop Classic	sonicaid	REPLACEMENT Newman Medical DigiDop II
NON DISPLAY D920 D930 Sonicaid D920/D930	DD-300	or	DD-330
SR2/SR3 FD1/FD3 Sonicaid SR2/SR3 FD1/FD3	DD-700	or	DD-770
DISPLAY SRX	DD-700	or	DD-770
SRXR Sonicaid SRX/SRXR	DD-701	or	DD-770R

PROBE ORDERING INFORMATION ON BACK SIDE




Newman Medical Replacement Probe Cross Reference

Vascular Probes



Fetal Probes



Professional and Sensitive Dopplers

Quality

- Better sound - Digital Sound Clarity™ produces up to 6x signal of other Dopplers
- Better warranty - 5 year honest warranty all parts & labor. Includes rechargeable batteries on rechargeable units.
- Better support & service - knowledgeable staff and guaranteed service

Sustainability

- Smart Recharge System™ (SRS)
 - Significant cost savings
 - New rechargeable batteries hold charge for years
 - May be recharged more than 2000 times

Innovation

- Continual product development and research

DigiDop Classic



- Standard size
- Longer battery life than similar Natus products

DigiDop II



- Smaller size
- Ergonomic - fits comfortably in your hand
- Similar battery life to Natus

DigiDop 901



- Ideal for hospital use
- Great room-filling sound from two speakers on top
- Two probe cradles
- Standard drop-on charging pyramid
- Optional wall mount recharger

How to build your Doppler

DOPPLER CODE

PROBE CODE

Doppler

- DD-300 Basic
- DD-301 Basic + Recharge
- DD-700 HR Display
- DD-701 HR Display + Recharge
- DD-901 Tabletop + HR Display + Recharge

Probe

- D2 2MHz fetal probe
- D2W 2MHz waterproof fetal probe
- D3 3MHz fetal probe
- D3W 3MHz waterproof fetal probe

DigiDop "Classic"

DigiDop II

- DD-330 Basic
- DD-330R Basic + Recharge
- DD-770 HR Display
- DD-770R HR Display + Recharge

- D5 5MHz vascular probe
- D8 8MHz vascular probe
- DPG Audio PPG probe - includes toe cuff

Interchangeable probes can be used on any DigiDop

Vascular



- Great sensitivity
- Ergonomic shape easy to hold, place in sterile glove finger for use in sterile vessels
- D8 - 8 MHz** - best for shallow vessels like radial, ankle, foot, digits, penile
- D5 - 5 MHz** - wider beam, deeper penetration. Good for brachial, carotid, femoral, deep veins

Fetal



- Great sensitivity
- Smooth shape, comfortable for patient, easy to hold and clean
- D3 - 3 MHz** - wide angle, ideal for early fetal hearts
- D3W** - a unique 3 MHz waterproof probe
- D2 - 2MHz** - deeper penetration
- D2W** - 2 MHz - waterproof

PPG



- DDPG**
- Unique audio PPG probe for use on digits. Great for toe-brachial studies

Arjo-Huntleigh-Newman Medical Part Number Cross-Reference Scan to learn more about DigiDop Doppler Systems

VASCULAR MODELS NUMBERS

Huntleigh	DD Classic	DD II	Huntleigh	DD Classic	DD II
D900-P-USA/VP8	DD-300-D8	DD-330-D8	LUS-DMX-VP10XS	DD-700-D8	DD-770-D8
LUS-D900-EZ8XS	DD-300-D8	DD-330-D8	LUS-DMX-VP4XS	DD-700-D5	DD-770-D5
LUS-D900-P-VP10KIT	DD-300-D8	DD-330-D8	LUS-DMX-VP5XS	DD-700-D5	DD-770-D5
LUS-D900-P-VP4KIT	DD-300-D5	DD-330-D5	LUS-DMX-VP8XS	DD-700-D8	DD-770-D8
LUS-D900-P-VP5KIT	DD-300-D5	DD-330-D5	LUS-SD2-EZ8XS	DD-700-D8	DD-770-D8
LUS-D900-VP10XS	DD-300-D5	DD-330-D5	LUS-SD2-P-EZ8KIT	DD-700-D8	DD-770-D8
LUS-D900-VP4XS	DD-300-D5	DD-330-D5	LUS-SD2-P-VP10KIT	DD-700-D8	DD-770-D8
LUS-D900-VP5XS	DD-300-D5	DD-330-D5	LUS-SD2-P-VP4KIT	DD-700-D5	DD-770-D5
LUS-D900-VP8XS	DD-300-D5	DD-330-D5	LUS-SD2-P-VP5KIT	DD-700-D5	DD-770-D5
LUS-DMX-EZ8XS	DD-700-D8	DD-770-D8	LUS-SD2-P-VP8KIT	DD-700-D8	DD-770-D8
LUS-DMXR-EZ8XS	DD-701-D8	DD-770R-D8	LUS-SD2-VP10XS	DD-700-D8	DD-770-D8
LUS-DMXR-VP10XS	DD-701-D8	DD-770R-D8	LUS-SD2-VP4XS	DD-700-D5	DD-770-D5
LUS-DMXR-VP4XS	DD-701-D5	DD-770R-D5	LUS-SD2-VP5XS	DD-700-D5	DD-770-D5
LUS-DMXR-VP5XS	DD-701-D5	DD-770R-D5	LUS-SD2-VP8XS	DD-700-D8	DD-770-D8
LUS-DMXR-VP8XS	DD-701-D8	DD-770R-D8			

FETAL MODEL NUMBERS

Huntleigh	DD Classic	DD II	Huntleigh	DD Classic	DD II
D920-P-USA	DD-300-D2W	DD-330-D2W	LUS-SRX-OP2XS	DD-700-D2	DD-770-D2
D920-USA	DD-300-D2W	DD-330-D2W	LUS-SRX-OP3XS	DD-700-D3	DD-770-D3
D930-P-USA	DD-300-D3W	DD-330-D3W	LUS-SRXR-OP2XS	DD-701-D2	DD-770R-D2
D930-USA	DD-300-D3W	DD-330-D3W	LUS-SRXR-OP3XS	DD-701-D3	DD-770R-D3
FD1-P-USA	DD-700-D2	DD-770-D2	SR2-USA	DD-700-D2W	DD-770-D2W
FD1-USA	DD-700-D2W	DD-770-D2W	SR3-USA	DD-700-D3W	DD-770-D3W
FD3-P-USA	DD-700-D3	DD-770-D3	Doppler Unit Only - No Probe		
FD3-USA	DD-700-D3W	DD-770-D3W	SRXR-USA	DD-701	DD-770R
			SRX-USA	DD-700	DD-770
			SD2-P-USA	DD-700	DD-770

