



Primary Care



Podiatry



Wound Care



Community Health

CLINIC PLACEMENT



Rapid ABI-Q

EXAMS PERFORMED



American-Made Quality



Includes labor, parts, and accidental breakage



ABI-Q Rapid Screen



DICOM Compatible

ABI-Q

The ABI-Q Rapid-Test System delivers fast, accurate PAD diagnosis in seconds, right in your office. Get comprehensive reports quickly, enabling early detection and improved patient outcomes. Cost-effective and portable, ABI-Q is ideal for primary care settings, especially for patients with diabetes or those with calcified arteries. Simplify testing, enhance care, and benefit from reimbursement opportunities with ABI-Q.

ABI-Q Exams

- ABI-Q Rapid Screen

ABI-Q Includes:

- Cuffs
 - 2 x 10cm EZ-Clean Vasocuffs
- Two-Year Warranty
 - Includes all parts, labor, and accidental breakage



simpleABI Vascular Systems: Quality You Can Trust

Made in the USA

- High-quality components, manufactured in Colorado.
- Durable design for busy clinical settings.
- 2-year warranty with complete coverage.



Newman Medical: Your Partner in Vascular Care

- Leaders in ultrasound diagnostics since the 1970s.
- Committed to innovation and providing exceptional vascular products.
- Focused on service and long-term customer relationships.

First-Class Support

- Experienced customer service team with extensive industry knowledge.
- Dedicated field representatives for personalized assistance.
- Comprehensive clinical support resources, both online and in-person.

ABI-Q



Exams Performed



ABI-Q Rapid Screen



Single-Level PVR

Ordering Information

System

SKU	DESCRIPTION
ABI-Q	simpleABI-Q Rapid Screen System

Accessories

SKU	DESCRIPTION
ABI-120	Microsoft® Windows® PC for use with simpleABI Systems
STND-130	simpleABI Roll Stand for Cuff-Link Systems
ACC-200	Set of 8 Cuff-Link Tubing
ACC-210	Set of 4 Cuff-Link Tubing
ACC-230	Memory Stick, USB with simpleABI Software
PUMP-100	Manual cuff inflator



Replacement Cuffs

SKU	DESCRIPTION
CUFF-200	EZ Clean Vascular Cuff, 10cm (ankle, small arm)
CUFF-210	EZ Clean Vascular Cuff, 12cm (arm, ankle, calf)
CUFF-210L	EZ Clean Vascular Cuff, 12cm long (above knee)
CUFF-220D	Digit Cuff, 2.5cm - Latex-Free
CUFF-240D	Digit Cuff, 1.9cm - Latex-Free
CUFF-SET-242	EZ Clean Vascular Cuff Set - 2x10cm, 4x12cm, 2x12cm - 8/pk

ABI-Q Technical Specifications

Weight:	3.6 kgs (8 pounds)
Level of Protection against electrical shock:	Type B Applied Part, Class II Equipment
Comply with the following Standards:	IEC60101-1, IEC60601-2
Transport/Storage Temperature:	-20° ~ 50°C (-4° ~ 122°F)
Transport/Storage Humidity:	5% ~ 90%, non-condensing
GTIN:	00852785008661

REIMBURSEMENT INFORMATION

Newman Medical's simpleABI vascular systems frequently utilize CPT codes 93922, 93923, and 93924 for billing and reimbursement purposes.

To obtain expert guidance on the appropriateness of these codes for specific systems, we highly recommend consulting with one of our experienced vascular specialists.

