

DigiDop

MADE IN DENVER, CO, USA

INNOVATIVE DOPPLERS FOR TODAY'S SAVVY PHYSICIAN

- Designed, manufactured & serviced in the U.S.
- Market leading 5-year warranty
- Rechargeable systems
- Extended battery life
- Interchangeable probes



THE NEWMAN MEDICAL DIFFERENCE

QUALITY

- Unprecedented 5-year warranty
 - includes labor, parts, breakage, and probes
- All probes tested
- Acoustical output power measured for all probes
- Certified to meet international quality standard ISO 13485

SUSTAINABILITY

- Rechargeable systems include premium batteries that can be charged 2,100+ times
- Smart Recharge System charges batteries fast and won't charge if alkaline batteries are accidentally inserted

INNOVATION

- Continual product development
- Investment in doppler innovation
- Advanced features include:
 - Digital Signal Clarity
 - Wide Angle XD 3 MHz probes
 - Smart Recharge System

HISTORY

- All Dopplers are designed, manufactured, and serviced in Colorado, U.S.
- The Newman family has been involved with Doppler since it was first introduced in the 1970s



MADE IN
DENVER, CO
UNITED STATES OF AMERICA

MADE IN

USA

INNOVATIVE TECHNOLOGY

Easy-to-Use | Dependable



THE DIGIDOP FAMILY OF PRODUCTS INCLUDES COUNTLESS TOP-OF-THE-LINE FEATURES:

5-year warranty covering parts, labor, and breakage
No other Doppler manufacturer comes close!



Large, bright display featuring battery condition, signal strength of Doppler or PPG signal, and volume setting (some models include calculation and display of autocorrelated heart rate)

Unbreakable, easy-to-clean magnetic cradle holds probes



Stainless steel microscreens protect speakers and keep dirt out

Premium batteries



Digital Signal Clarity sound processing

Interchangeable probes for vascular and obstetrical use

NEWMAN MEDICAL DOPPLERS



DIGIDOP

- Standard size
- Superior battery endurance compared to competitors
- Robust & solid
- Large, bright display



DD 300|301|700|701

DIGIDOP II

- Compact size
- Comfortable palm-fitting design
- Equivalent battery endurance to competitors
- Large, bright display



DD 330|330R|770|770R

DIGIDOP 901

- Ideal for hospital use
- Clear stereo sound from dual top speakers
- Two probe cradles
- Standard drop-on charging pyramid
- Optional wall mount recharger



DD 901

How to build your Doppler | Doppler

Doppler Code – Probe Code

Example

DD-901-D2

| Probe

DD300	Basic	D2	2MHz obstetrical probe
DD301	Rechargeable	D2W	2MHz waterproof obstetrical probe
DD700	Includes heartrate	D3	3MHz obstetrical probe
DD701	Rechargeable, HR	D3W	3MHz waterproof obstetrical probe
DD330	Basic	D5	5MHz vascular probe
DD330R	Rechargeable	D8	8MHz vascular probe
DD770	Includes heartrate	DPPG	Audio PPG Probe - Includes toe cuff
DD770R	Rechargeable, HR		
DD901	Rechargeable, HR		

INTERCHANGEABLE PROBES CAN BE USED ON ANY DIGIDOP



VASCULAR

- Great sensitivity
- Ergonomic shape easy to hold, place in sterile glove finger for use in sterile fields
- **D8 – 8 MHz**
Best for shallow vessels like radial, ankle, foot, digits, penile
- **D5 – 5 MHz**
Wider beam, deeper penetration. Good for brachial, carotid, femoral, deep veins



FETAL

- Great sensitivity
- Smooth shape comfortable for patient, easy to hold and clean
- **D3 – 3 MHz** - wide angle, ideal for early fetal hearts
- **D3W** - a unique 3 MHz waterproof probe
- **D2 – 2MHz**
deeper penetration
- **D2W – 2 MHz**
waterproof



PPG

- **DPPG** - Unique audio PPG probe for use on digits
- Great for toe-brachial studies

800-267-5549

info@newman-medical.com
newman-medical.com

Newman  Medical

Employee Owned & Operated

POPULAR ACCESSORIES

For a complete listing of accessories go to newman-medical.com or call customer service

MADE IN

USA



WALL MOUNT



ROLL STANDS



ABI-TBI KITS

Newman  Medical

MADE IN
DENVER, CO
UNITED STATES OF AMERICA



800-267-5549

info@newman-medical.com
newman-medical.com



DROP-IN
HOLDER/
CHARGER/

Smart  Recharge
System