


Newman Medical

Product Catalog

MADE IN
USA

simple ABI



DigiDop



Automated Cuff-Link Systems

Advanced Automated Testing with Dual Capabilities

Our **Cuff-Link systems** combine automated efficiency with comprehensive testing capabilities. All models feature **dual testing**: rapid ABI-Q screening AND comprehensive Doppler-based exams, with automated cuff selection and inflation to reduce user error.

Key Features

- Dual testing capabilities (ABI-Q + Doppler)
- Automated cuff selection & inflation
- Streamlined workflow
- Reduced user error
- Variable testing options by model
- Adjustable-height roll stand included
- Advanced diagnostic capabilities

2-YEAR INDUSTRY-LEADING WARRANTY

400CL
ABI + TBI
4-Port System

450CL
+ Stress/Exercise
4-Port System

500CL
+ Segmental
8-Port System

600CL
Complete Testing
8-Port System



simpleABI-Q Rapid-Test System

Doppler-Free PVR-Based Rapid Testing

The **simpleABI-Q** is designed for busy practices needing fast, efficient PAD screening with limited time, staff, or resources. This innovative system uses PVR waveforms and proprietary software algorithms to provide accurate PAD assessments in just **1-2 minutes** without requiring Doppler skills.

Key Features

- Fast 1-2 minute exams
- Push-button operation
- No Doppler skills required
- Perfect for diabetics with calcified arteries
- Portable design
- Minimal training required
- Cost-effective screening
- 2-year industry-leading warranty

2-YEAR INDUSTRY-LEADING WARRANTY



Available Exams:
Rapid ABI-Q | ABI-Q + Exercise
2-Year Industry-Leading Warranty

Manual Handheld Doppler-Based Systems

Budget-Friendly Manual Cuff Inflation with Doppler Technology

Our **manual systems** provide reliable, accurate ABI measurements using high-quality Doppler technology with manual cuff inflation. Perfect for budget-conscious practices seeking foundational to comprehensive PAD testing capabilities

Key Features

- Portable and compact design
- Manual cuff inflation control
- High-quality Doppler technology
- Cost-effective solution
- Made in USA
- Minimal training required
- 2-year industry-leading warranty
- Optional SteadyDop probe available

2-YEAR INDUSTRY-LEADING WARRANTY



simpleABI-250
Basic PAD Screening
Doppler ABI

simpleABI-300
ABI + TBI Testing
DPPG Probe Included



Professional & Sensitive Doppler Solutions
Interchangeable Fetal & Vascular Probes

MADE IN
USA

DigiDop Classic

Standard Professional Size | Extended Battery Life | Superior Performance

Key Features

- Standard professional size for optimal handling
- Digital Sound Clarity™ technology
- Extended battery life
- Comfortable grip for extended sessions
- Robust construction for daily clinical use

Available Options

- Non-Display: DD-300 (Alkaline) or DD-301 (Rechargeable)
- Heart Rate Display: DD-700 (Alkaline) or DD-701 (Rechargeable)
- Optional recharge base for rechargeable models

Available Models: DD-300 | DD-301 | DD-770 | DD-700 | DD-701



Quickly switch between fetal and vascular applications with interchangeable probes.

DigiDop II

Ultra-Compact Pocket-Sized | Maximum Portability | Mobile Healthcare

Key Features

- Ultra-compact, pocket-sized design
- Ergonomic shape for comfortable handling
- Maximum portability without sacrificing performance
- Digital Sound Clarity™ technology
- Perfect for mobile healthcare professionals

Available Options

- Non-Display: DD-330 (Alkaline) or DD-330R (Rechargeable)
- Heart Rate Display: DD-770 (Alkaline) or DD-770R (Rechargeable)
- Audio Recording: DD-330AR (Rechargeable) or DD-770AR available (Rechargeable)
- Optional recharge base for rechargeable models

Available Models: DD-330 | DD-330R | DD-770 | DD-770R | DD-330AR | DD-770AR



Quickly switch between fetal and vascular applications with interchangeable probes.

DigiDop Tabletop (DD-901)

Hospital Grade | Room-Filling Sound | Dual Probe Cradles

Key Features

- Dual speakers for superior sound projection
- Two convenient probe cradles
- Drop-on charging pyramid
- Optional wall mount charging capability
- Ideal for hospital and clinical environments
- Heart rate display and rechargeable

Clinical Benefits

- Room-filling stereo sound quality
- Hands-free operation with cradles
- Always charged and ready to use
- Professional stationary solution
- Gold standard in clinical diagnostics



Quickly switch between fetal and vascular applications with interchangeable probes.

5-YEAR INDUSTRY-LEADING WARRANTY • ALL PARTS & LABOR • NO HIDDEN FEES • MADE IN USA

800-267-5549 | info@newman-medical.com | www.newman-medical.com

PROBE SOLUTIONS - Precision diagnostic probes for every clinical application

simpleABI[®] System Probes

NEW SteadyDop Technology
Revolutionary stethoscope-like grip eliminates movement errors with wide sensors for easy artery tracking. Built-in optimal angles ensure perfect positioning for instant detection every time.



8MHz General Purpose Probe



D8

Optimized for peripheral vessels in vascular applications. Our most versatile probe for routine ABI testing.

STANDARD Included with all systems

Compatible: All simpleABI manual & automated systems

5MHz Deep Vessel Probe



D5

Optimized for deeper vessels in challenging vascular applications where standard probes may struggle.

OPTIONAL Enhanced capability

Compatible: All simpleABI manual & automated systems

8MHz SteadyDop Probe



D8

Wide sensors with revolutionary grip design for peripheral vessels. Eliminates probe movement errors for consistent readings.

NEW PREMIUM OPTIONAL

Compatible: All simpleABI manual & automated systems

5MHz SteadyDop Probe



D5

Wide sensors with revolutionary grip design for deeper vessels. Perfect for challenging anatomies and extended assessments.

NEW PREMIUM OPTIONAL

Compatible: All simpleABI manual & automated systems

DPPG Audio Probe



DPPG

For incompressible arteries when standard Doppler is insufficient. Provides both audio and visual PPG waveforms.

SPECIALIZED

Compatible: Cuff-Link 3500+ series & ABI-300 manual

SPPG Probe



SPPG

Connects directly to Cuff-Link control unit for PPG waveforms during TBI exams. Advanced diagnostic capability.

TBI SPECIALIST

Compatible: Cuff-Link 4000+ series only

DigiDop[®] System Probes

3MHz General Purpose Fetal



D3

Optimized for early fetal heartbeat detection and use throughout pregnancy. Most versatile fetal probe.

STANDARD

Compatible: All DigiDop models

2MHz Later-Term Fetal



D2

Optimized for later term pregnancies or larger patients where deeper penetration is needed.

SPECIALIZED

Compatible: All DigiDop models

3MHz Waterproof Fetal



D3W

Designed for water applications and labor deliveries. Fully submersible for water birth monitoring.

WATER BIRTH

Compatible: All DigiDop models

2MHz Waterproof Fetal



D2W

Waterproof version for larger patients during water applications and labor deliveries.

WATER BIRTH

Compatible: All DigiDop models

8MHz Peripheral Vascular



D8

General purpose vascular probe optimized for peripheral vessels. Ideal for routine vascular assessments.

STANDARD Included with all systems

Compatible: All DigiDop models

5MHz Deep Vascular



D5

Optimized for deeper vessels in vascular applications where standard frequency is insufficient.

OPTIONAL Enhanced capability

Compatible: All DigiDop models

8MHz SteadyDop Probe



D8

Wide sensors with stethoscope-like grip for peripheral vessels. Enhanced stability and accuracy.

NEW PREMIUM OPTIONAL

Compatible: All DigiDop models

5MHz SteadyDop Probe



D5

Wide sensors with revolutionary grip for deeper vessels. Perfect for challenging vascular assessments.

NEW PREMIUM OPTIONAL

Compatible: All DigiDop models

DPPG Audio Probe



DPPG

For incompressible arteries when standard Doppler probes are insufficient. Specialized diagnostic capability.

SPECIALIZED

Compatible: All DigiDop models

Newman Medical

www.newman-medical.com
800-267-5549 | info@newman-medical.com

FETAL

VASCULAR