

Edan Warranty

- 2 Year warranty on the Sonotrax Doppler main unit
- 1 Year warranty on the accessories or disposables

DigiDop

EDANUSA

Model

Newman Medical
DigiDop Classic

Newman Medical
DigiDop II



**Sonotrax
Lite**

DD-300



or

DD-330



**Sonotrax
Vascular**

DD-300



or

DD-330



**Sonotrax
Basic A**

DD-700



or

DD-770



**Sonotrax
Pro**

N/A

or

DD-770A



**Sonotrax
II**

DD-701
⚡



or

DD-770R
⚡



**Sonotrax
II Pro**

N/A

or

DD-330AR
⚡



EDANUSA

Model



**SD5/SD6
Ultrasonic
Tabletop
Doppler**



DD-901 ⚡

NON-DISPLAY

HANDHELD DOPPLERS

DISPLAY

TABLETOP



Includes labor, parts,
and accidental breakage



American-Made
Quality



Includes cords
and probes



Waterproof
Option



Rechargeable
Option

PROBE ORDERING INFORMATION ON BACK SIDE



Newman Medical Replacement Probe Cross Reference

Vascular Probes



Model D8
8MHz Vascular Probe

or

Model D5
5MHz Vascular Probe



Model D8-SD
8 MHz SteadyDop Vascular Probe

or

Model D5-SD
5 MHz SteadyDop Vascular Probe



Model DPPG
Audio TBI PPG Probe



Model D2
2MHz Fetal Probe

or

Model D3
3MHz Fetal Probe

Fetal Probes

Model D2W
2MHz Waterproof Fetal Probe

or

Model D3W
3MHz Waterproof Fetal Probe



Professional and Sensitive Dopplers

Quality

- Better sound - Digital Sound Clarity™ produces up to 6x signal of other Dopplers
- Better warranty - 5 year honest warranty all parts & labor. Includes rechargeable batteries on rechargeable units.
- Better support & service - knowledgeable staff and guaranteed service

Sustainability

- Smart Recharge System™ (SRS)
 - Significant cost savings
 - New rechargeable batteries hold charge for years
 - May be recharged more than 2000 times

Innovation

- Continual product development and research

DigiDop Classic



- Standard size
- Longer battery life than similar Natus products

DigiDop II



- Smaller size
- Ergonomic - fits comfortably in your hand
- Similar battery life to Natus

DigiDop 901



- Ideal for hospital use
- Great room-filling sound from two speakers on top
- Two probe cradles
- Standard drop-on charging pyramid
- Optional wall mount recharger

How to build your Doppler

DOPPLER CODE

PROBE CODE

Doppler

DD-300	Basic
DD-301	Basic + Recharge
DD-700	HR Display
DD-701	HR Display + Recharge
DD-901	Tabletop + HR Display + Recharge

DigiDop "Classic"

DD-330	Basic
DD-330R	Basic + Recharge
DD-770	HR Display
DD-770R	HR Display + Recharge

DigiDop II

Probe

D2	2MHz fetal probe
D2W	2MHz waterproof fetal probe
D3	3MHz fetal probe
D3W	3MHz waterproof fetal probe
D5	5MHz vascular probe
D8	8MHz vascular probe
DPPG	Audio PPG probe - includes toe cuff

Interchangeable probes can be used on any DigiDop

Vascular



- Great sensitivity
- Ergonomic shape easy to hold, place in sterile glove finger for use in sterile fields
- D8 - 8 MHz** - best for shallow vessels like radial, ankle, foot, digits, penile
- D5 - 5 MHz** - wider beam, deeper penetration. Good for brachial, carotid, femoral, deep veins

Fetal



- Great sensitivity
- Smooth shape, comfortable for patient, easy to hold and clean
- D3 - 3 MHz** - wide angle, ideal for early fetal hearts
- D3W** - a unique 3 MHz waterproof probe
- D2 - 2MHz** - deeper penetration
- D2W** - 2 MHz - waterproof

PPG



- DPPG**
- Unique audio PPG probe for use on digits. Great for toe-brachial studies

Edan Diagnostics - DigiDop Doppler SKU Cross-Reference

SonoTrax Vascular

Edan	DigiDop Classic	DigiDop II
STV.04	DD-300-D5	DD-330-D5
STV.05	DD-300-D5	DD-330-D5

SonoTrax Lite Fetal

Edan	DigiDop Classic	DigiDop II
STL.02	DD-300-D2	DD-330-D2
STL.03	DD-300-D3	DD-330-D3

SonoTrax Basic A

Edan	DigiDop Classic	DigiDop II
STBA.02	DD-700-D2	DD-770-D2
STBA.03	DD-700-D3	DD-770-D3
STBA.04	DD-700-D5	DD-770-D5
STBA.05	DD-700-D5	DD-770-D5
STBA.08	DD-700-D8	DD-770-D8

SonoTrax II

Edan	DigiDop Classic	DigiDop II
ST2.02	DD-701-D3	DD-770R-D3
ST2.03	DD-701-D5	DD-770R-D5
ST2.04	DD-701-D5	DD-770R-D5
ST2.05	DD-701-D5	DD-770R-D5
ST2.08	DD-701-D8	DD-770R-D8

SonoTrax Pro

Edan	DigiDop Classic	DigiDop II
STP.02	N/A	DD-770A-D2
STP.03	N/A	DD-770A-D3
STP.04	N/A	DD-770A-D5
STP.05	N/A	DD-770A-D5
STP.08	N/A	DD-770A-D8

SonoTrax Pro II

Edan	DigiDop Classic	DigiDop II
ST2P.02	N/A	DD-330AR-D2
ST2P.03	N/A	DD-330AR-D3
ST2P.04	N/A	DD-330AR-D5
ST2P.05	N/A	DD-330AR-D5
ST2P.08	N/A	DD-330AR-D8

Newman Medical



MADE IN
USA

Scan to learn more
about DigiDop Dopplers



Visit www.newman-medical.com or call 800-267-5549 for more information