



MADE IN DENVER, CO, USA

INNOVATIVE DOPPLERS FOR TODAY'S SAVVY PHYSICIAN

- Designed, manufactured & serviced in the U.S.
- Market leading 5-year warranty
- Rechargeable systems
- Extended battery life
- Interchangeable probes



Newman  Medical

Employee Owned & Operated

THE NEWMAN MEDICAL DIFFERENCE

QUALITY

- Unprecedented 5-year warranty
 - includes labor, parts, breakage, and probes
- All probes tested
- Acoustical output power measured for all probes
- Certified to meet international quality standard ISO 13485

SUSTAINABILITY

- Rechargeable systems include premium batteries that can be charged 2,100+ times
- Smart Recharge System charges batteries fast and won't charge if alkaline batteries are accidentally inserted

INNOVATION

- Continual product development
- Investment in doppler innovation • Advanced features include:
 - Digital Signal Clarity,
 - Wide Angle XD 3 MHz probes,
 - Smart Recharge System

HISTORY

- All Dopplers are designed, manufactured, and serviced in Colorado, U.S.
- The Newman family has been involved with Doppler since it was first introduced in the 1970s



MADE IN
DENVER, CO
UNITED STATES OF AMERICA

INNOVATIVE TECHNOLOGY

Easy-to-Use | Dependable



THE DIGIDOP FAMILY OF PRODUCTS INCLUDES COUNTLESS TOP-OF-THE-LINE FEATURES:

5-year warranty covering parts, labor, and breakage
No other Doppler manufacturer comes close!

Large, bright display featuring battery condition, signal strength of Doppler or PPG signal, and **volume setting** (some models include calculation and display of autocorrelated heart rate)



Unbreakable, easy-to-clean magnetic cradle holds probes



Stainless steel microscreens protect speakers and keep dirt out

Premium batteries



Digital Signal Clarity
sound processing

Interchangeable probes for
vascular and obstetrical use

NEWMAN MEDICAL DOPPLERS



DIGIDOP

- Standard size
- Superior battery endurance compared to competitors
- Robust & solid
- Large, bright display



DD 300|301|700|701



DIGIDOP II

- Compact size
- Comfortable palm-fitting design
- Equivalent battery endurance to competitors
- Large, bright display



DD 330|330R|770|770R



DIGIDOP 901

- Ideal for hospital use
- Clear stereo sound from dual top speakers
- Two probe cradles
- Standard drop-on charging pyramid
- Optional wall mount recharger

DD 901

How to build your Doppler | Doppler

Doppler Code – Probe Code

Example

DD-901-D2

Probe

DD300	Basic
DD301	Rechargeable
DD700	Includes heartrate
DD701	Rechargeable, HR
DD330	Basic
DD330R	Rechargeable
DD770	Includes heartrate
DD770R	Rechargeable, HR
DD901	Rechargeable, HR

D2	2MHz obstetrical probe
D2W	2MHz waterproof obstetrical probe
D3	3MHz obstetrical probe
D3W	3MHz waterproof obstetrical probe
D5	5MHz vascular probe
D8	8MHz vascular probe
DPPG	Audio PPG Probe - Includes toe cuff

INTERCHANGEABLE PROBES CAN BE USED ON ANY DIGIDOP

MADE IN

USA



VASCULAR

- Great sensitivity
- Ergonomic shape easy to hold, place in sterile glove finger for use in sterile fields
- **D8 – 8 MHz**
Best for shallow vessels like radial, ankle, foot, digits, penile
- **D5 – 5 MHz**
Wider beam, deeper penetration. Good for brachial, carotid, femoral, deep veins



OBSTETRICAL

- Great sensitivity
- Smooth shape comfortable for patient, easy to hold and clean
- **D3 – 3 MHz** - wide angle, ideal for early fetal hearts
- **D3W** - a unique 3 MHz waterproof probe
- **D2 – 2MHz**
deeper penetration
- **D2W – 2 MHz**
waterproof



PPG

- **DPPG** - Unique audio PPG probe for use on digits
- Great for toe-brachial studies

800-267-5549

info@newman-medical.com
newman-medical.com

Newman  Medical

Employee Owned & Operated

POPULAR ACCESSORIES

For a complete listing of accessories go to newman-medical.com or call customer service



WALL MOUNT



ROLL STANDS



ABI-TBI KITS

Newman  Medical

Employee Owned & Operated

MADE IN
DENVER, CO
UNITED STATES OF AMERICA



800-267-5549

info@newman-medical.com
newman-medical.com



DROP-IN
HOLDER/
CHARGER/