

MADE IN DENVER, CO, USA

# INNOVATIVE DOPPLERS FOR TODAY'S SAVVY PHYSICIAN

- Designed, manufactured & serviced in the U.S.
- Market leading 5-year warranty
- Rechargeable systems
- Extended battery life
- Interchangeable probes







### THE **NEWMAN MEDICAL** DIFFERENCE

#### **QUALITY**

- Unprecedented 5-year warranty
  - o includes labor, parts, breakage, and probes
- All probes tested
- Acoustical output power measured for all probes
- Certified to meet international quality standard ISO 13485

#### **SUSTAINABILITY**

- Rechargeable systems include premium batteries that can be charged 2,100+ times
- Smart Recharge System charges batteries fast and won't charge if alkaline batteries are accidentally inserted

#### **INNOVATION**

- Continual product development
- Investment in doppler innovation
- Advanced features include:
  - Digital Signal Clarity
  - Wide Angle XD 3 MHz probes
  - Smart Recharge System

#### **HISTORY**

- All Dopplers are designed, manufactured, and serviced in Colorado, U.S.
- The Newman family has been involved with Doppler since it was first introduced in the 1970s

3 XD **MADE IN** 

DENVER, CO
UNITED STATES OF AMERICA

# INNOVATIVE TECHNOLOGY Easy-to-Use | Dependable



THE DIGIDOP FAMILY OF PRODUCTS INCLUDES COUNTLESS TOP-OF-THE-LINE FEATURES:

## **5-year warranty covering parts, labor, and breakage**No other Doppler manufacturer comes close!



Large, bright display featuring battery condition, signal strength of Doppler or PPG signal, and volume setting (some models include calculation and display of autocorrelated heart rate)



Stainless steel microscreens protect speakers and keep dirt out

Digital Signal Clarity sound processing

Unbreakable, easy-to-clean magnetic cradle holds probes



**Premium** batteries



Interchangeable probes for vascular and obstetrical use

### NEWMAN MEDICAL DOPPLERS





#### **DIGIDOP**

- Standard size
- Superior battery endurance compared to competitors
- Robust & solid
- Large, bright display



DD 300|301|700|701



#### **DIGIDOP II**

- Compact size
- Comfortable palm-fitting design
- Equivalent battery endurance to competitors
- Large, bright display



DD 330|330R|770|770R



#### DIGIDOP 901

- Ideal for hospital use
- Clear stereo sound from dual top speakers
- Two probe cradles
- Standard drop-on charging pyramid
- Optional wall mount recharger



DD 901

#### How to build your Doppler | Doppler

Doppler Code – Probe Code

Example DD-901-D2

DD300	Basic
DD301	Rechargeable
DD700	Includes heartrate
DD701	Rechargeable, HR
DD330	Basic
DD330R	Rechargeable
DD770	Includes heartrate
DD770R	Rechargeable, HR
DD901	Rechargeable HR

#### **Probe**

D2	2MHz obstetrical probe
D2W	2MHz waterproof obstetrical probe
D3	3MHz obstetrical probe
D3W	3MHz waterproof obstetrical probe
D5	5MHz vascular probe
D8	8MHz vascular probe
DPPG	Audio PPG Probe - Includes toe cuff

# INTERCHANGEABLE PROBES CAN BE USED ON ANY DIGIDOP





#### **VASCULAR**

- Great sensitivity
- Ergonomic shape easy to hold, place in sterile glove finger for use in sterile fields
- D8 8 MHz

Best for shallow vessels like radial, ankle, foot, digits, penile

• D5 - 5 MHz

Wider beam, deeper penetration. Good for brachial, carotid, femoral, deep veins



#### **FETAL**

- Great sensitivity
- Smooth shape comfortable for patient, easy to hold and clean
- D3 3 MHz wide angle, ideal for early fetal hearts
- **D3W** a unique 3 MHz waterproof probe
- **D2** 2MHz deeper penetration
- **D2W** 2 MHz waterproof



#### **PPG**

- **DPPG** Unique audio PPG probe for use on digits
- Great for toe-brachial studies

800-267-5549

info@newman-medical.com newman-medical.com



### POPULAR ACCESSORIES

For a complete listing of accessories go to newmanmedical.com or call customer service





**WALL MOUNT** 



**ROLL STANDS** 



**ABI-TBI KITS** 



DENVER, CO
UNITED STATES OF AMERICA



800-267-5549

info@newman-medical.com newman-medical.com

