

DIGIDOP DD-901 TABLETOP

Heart Rate Display, Rechargeable, Tabletop, Dual Probe Holder, Room-Filling Sound



VASCULAR



HR DISPLAY



RECHARGEABLE

Newman Medical



American-Made
Quality



Includes labor, parts,
and accidental
breakage



Includes cords
and probes



Rechargeable

Vascular | HR-Display Doppler System



The DigiDop Tabletop Doppler offers a feature-rich experience for healthcare professionals. This innovative system features Digital Sound Clarity™ technology for clear and accurate readings. Dual speakers provide room-filling sound for easy sharing of results. Two probe cradles ensure convenient probe storage, while a standard drop-on charging pyramid with a long battery life and an optional wall mount recharger offer increased flexibility.

The DigiDop Tabletop Doppler is built to last. Manufactured in the USA with high-quality materials, each DigiDop undergoes rigorous testing by Newman Medical to guarantee unwavering performance and durability. Invest in the DigiDop Tabletop Doppler and experience a dependable medical device with exceptional features.

DigiDop Doppler - The Clear Choice in Diagnostics

Superior Performance & Reliability:

- **Exceptional Clarity:** Digital Sound Clarity DSC™ technology enhances signals up to 6x, delivering superior diagnostic accuracy.
- **American Craftsmanship:** Proudly designed, manufactured, and serviced in the U.S., ensuring premium quality and rapid support.
- **Durability By Design:** Features an unbreakable probe holder built for years of dependable clinical use.
- **Versatile Functionality:** All probes are fully interchangeable, maximizing your investment and clinical capabilities.
- **Rechargeable System Option:** Eliminates ongoing battery expenses and reduces environmental impact.

Industry-Leading Warranty:

- **Comprehensive 5-Year Coverage:** Full parts and labor warranty with zero hidden fees or diagnostic charges.

Newman Medical: Pioneers in Ultrasound Excellence:

- **50+ Years of Innovation:** The Newman family has advanced Doppler technology since the 1970s.
- **Proven Industry Leadership:** Founders behind respected brands including Imex Medical Systems (Nicolet/Natus) and Summit Doppler.
- **Continuous Advancement:** Ongoing development of next-generation diagnostic technologies.
- **Partnership Approach:** Committed to your long-term success with exceptional service and support.

DIGIDOP DD-901 TABLETOP

VASCULAR | TABLETOP DOPPLER SYSTEM

Newman Medical

Accessories



Wall Board Mount
(ACC-251)



Roll Stand
(STND-150)



Ultrasound Gel
8.8oz | 2.1oz
(GEL-110) (GEL-100)



Ask about our Doppler Recycling Program

NEW



SIMPLIFY VASCULAR ASSESSMENTS WITH THE STEADYDOP PROBE

Its unique stethoscope-like grip and optimized design ensure easy vessel tracking. Compatible with Newman Medical's DigiDop & simpleABI systems. Optional accessory.

Contact your Newman Medical rep for details.

Ordering Information

Doppler System

SKU	DESCRIPTION
DD-901	Vascular Tabletop Doppler System

Interchangeable Probe Choices

SKU	DESCRIPTION
D5	5 MHz Vascular
D5-SD	5 MHz SteadyDop Vascular
D8	8 MHz Vascular
D8-SD	8 MHz SteadyDop Vascular
D2	2 MHz Fetal
D3	3 MHz Fetal



Accessories

SKU	DESCRIPTION
ACC-251	Wall Mount
STND-150	Roll stand with basket
GEL-110	Ultrasound gel 8.8oz 12/box
GEL-100	Ultrasound gel 2.1oz 12/box

How to order your Doppler

DOPPLER SYSTEM CODE (DD-###) - PROBE CODE (D#)

Example: DD-770-D3

DD-901 Technical Specifications

Dimensions:	210 x 155 x 115 mm (8.25 x 6 x 4.5 in)
Weight:	850 grams (30 ounces)
Battery Type:	3 x 1.5 nominal volts (AA/R6)
Level of Protection against electrical shock:	Type B Applied Part, Class II Equipment
Comply with the following Standards:	IEC60101-1, IEC60601-2, IEC60601-2-37
Operating Temperature:	10° ~ 40°C (50° ~ 104°F)
Operating Humidity:	30% ~ 75%
Transport/Storage Temperature:	-20° ~ 50°C (-4° ~ 122°F)
Transport/Storage Humidity:	5% ~ 90%, non-condensing
Battery Life (AA/R6):	600 mins (NiMH), 800 mins (Alkaline) 50 ~ 220
Heart Rate Range:	BPM
Heart Rate Calculation Accuracy:	±3 BPM
Audio Output Power	1.1W
GTIN:	852785008982

The DigiDop PRODUCT FAMILY

DigiDop "Classic"



- Standard size
- Longer battery life than similar competitor products

DigiDop II



- Smaller size
- Ergonomic - fits comfortably in your hand
- Similar battery life to competitor products

DigiDop 901 Tabletop

- Ideal for hospital use
- Great room-filling sound from two speakers on top
- Two probe cradles
- Standard drop-on charging pyramid
- Optional wall mount recharger

