

ENHANCE COMFORT AND MOBILITY

DigiDop Adjustable-Height Mobile Stand

The DigiDop Mobile Stand transforms your Doppler into a flexible, fully portable workstation, perfect for busy clinic environments, exam rooms, or shared workspaces. Designed with both the patient and operator in mind, this stand maximizes comfort, minimizes strain, and keeps your DigiDop system ready to go where you need it.



KEY FEATURES & CLINICAL BENEFITS

- **Adjustable Height for Comfort:** Easily adjust the stand height to accommodate patients of varying sizes and allow operators to work from a comfortable standing or sitting position, reducing back and neck strain.
- **Maximum Mobility:** Equipped with smooth-rolling casters (two locking, three non-locking) for quick, quiet movement between rooms. Lock the wheels securely for stable operation during exams.
- **Dedicated Storage Basket:** An integrated basket provides immediate, convenient storage for probes, gel, cables, and other supplies, keeping everything you need within reach.
- **Durable Construction:** Built with a sturdy aluminum alloy base and a stainless steel top pole, ensuring long-lasting durability and stability in demanding clinical settings.

The DigiDop PRODUCT FAMILY

DigiDop "Classic"

- Standard size
- Longer battery life than similar competitor products



DigiDop II

- Smaller size
- Ergonomic - fits comfortably in your hand
- Similar battery life to competitor products



DigiDop 901 Tabletop

- Ideal for hospital use
- Great room-filling sound from two speakers on top
- Two probe cradles
- Standard drop-on charging pyramid
- Optional wall mount recharger



Mobile Stand Specifications

Height adjustment	700~1020mm
Handle	Optional
Top pole	Φ25±0.5mm/Stainless steel
Lower pole	Φ50±1mm/Aluminum Alloy
Base	Φ520mm /Aluminum Alloy
Casters	two locking and three non-locking
Cart maximum load	7Kg
Basket	300*213*150mm / ABS
Basket maximum load	2kg
Desktop size	400*420mm

Ordering Information

DigiDop Doppler Model	Compatible Mobile Stand
DD-300	STND-101
DD-700	STND-101
DD-301	STND-111
DD-701	STND-111
DD-330	STND-121
DD-770	STND-121
DD-330R	STND-131
DD-770R	STND-131
DD-901	STND-150