

# DIGIDOP PROBES

FOR DIGIDOP CLASSIC AND DIGIDOP II DOPPLER SYSTEMS



American-Made  
Quality



Includes labor, parts,  
and accidental  
breakage



## D2

The DigiDop 2MHz fetal probe is optimized for later term pregnancies or larger patients.



## D2W

The DigiDop 2MHz waterproof fetal probe is designed for water applications and labor deliveries.



## D3

The DigiDop 3MHz probe is our general purpose obstetrical probe. It is optimized for early fetal heartbeat detection and use throughout pregnancy.



## D3W

The DigiDop 3MHz waterproof probe is designed for water applications and labor deliveries.



## D5

The DigiDop 5MHz probe is optimized for deeper vessels in vascular applications.



## D8

The DigiDop 8MHz probe is our general purpose vascular probe. The DigiDop 8MHz probe is optimized for peripheral vessels in vascular applications.



## DD-PPG

The DigiDop PPG probe is for use in vascular applications. The DigiDop PPG probe is for use when a Doppler probe is insufficient – i.e. for those with incompressible arteries.