




Primary Care



Podiatry



Wound Care



Single-Level



TBI w/PPG

CLINIC PLACEMENT

EXAMS PERFORMED



American-Made Quality



Includes labor, parts, and accidental breakage



DICOM Compatible

simpleABI

The simpleABI-300 Handheld PVR system is a portable and cost-effective solution for the accurate detection of Peripheral Artery Disease (PAD). This clinically precise system offers enhanced manual control, simplifies reimbursement, and is compact enough for diverse healthcare settings, including mobile services and community health fairs.

Proudly made in the USA with the highest quality standards, simpleABI is backed by a 2-year warranty for added peace of mind.

ABI-300 Exams

- Single-Level
- TBI with PPG

ABI-300 Includes:

- simpleABI software
- Doppler with 8 MHz probe
- Manual cuff inflator
- DPPG probe
- PVR
- Cuffs
 - 2 x 10cm EZ-Clean Vasocuffs
 - 2 x 12cm EZ-Clean Vasocuffs
 - 1 x 1.9cm digit cuff
- Carry bag
- User manual / Quick reference guide
- Two-Year Warranty
 - Includes all parts, labor, and accidental breakage



simpleABI Vascular Systems: Quality You Can Trust

Made in the USA

- High-quality components, manufactured in Colorado.
- Durable design for busy clinical settings.
- 2-year warranty with complete coverage.

Newman Medical: Your Partner in Vascular Care

- Leaders in ultrasound diagnostics since the 1970s.
- Committed to innovation and providing exceptional vascular products.
- Focused on service and long-term customer relationships.

First-Class Support


- Experienced customer service team with extensive industry knowledge.
- Dedicated field representatives for personalized assistance.
- Comprehensive clinical support resources, both online and in-person.



ABI-300



Exams Performed



Ordering Information

System

SKU	DESCRIPTION
ABI-300	simpleABI Single-Level Cuff-Link System with TBI

Accessories

SKU	DESCRIPTION
ABI-120	Microsoft® Windows® PC for use with simpleABI Systems
STND-130	simpleABI Roll Stand for Manual Systems
ACC-175	simpleABI-300 carry bag
PUMP-100	Manual cuff inflator

Replacement Cuffs

SKU	DESCRIPTION
CUFF-200	EZ Clean Vascular Cuff, 10cm (ankle, small arm)
CUFF-210	EZ Clean Vascular Cuff, 12cm (arm, ankle, calf)
CUFF-210L	EZ Clean Vascular Cuff, 12cm long (above knee)
CUFF-220D	Digit Cuff, 2.5cm - Latex-Free
CUFF-240D	Digit Cuff, 1.9cm - Latex-Free
CUFF-SET-242	EZ Clean Vascular Cuff Set - 2x10cm, 4x12cm, 2x12cm - 8/pk



ABI-300 Technical Specifications

Weight:	3.6 kgs (8 pounds)
Level of Protection against electrical shock:	Type B Applied Part, Class II Equipment
Comply with the following Standards:	IEC60101-1, IEC60601-2, IEC60601-2-37
Operating Temperature:	10° ~ 40°C (50° ~ 104°F)
Operating Humidity:	30% ~ 75%
Transport/Storage Temperature:	-20° ~ 50°C (-4° ~ 122°F)
Transport/Storage Humidity:	5% ~ 90%, non-condensing
Battery Voltage, Type (Doppler Unit):	3 x 1.5 volts, AA
Battery Life::	800, 1 minute exams
GTIN:	00852785008647

REIMBURSEMENT INFORMATION

Newman Medical's simpleABI vascular systems frequently utilize CPT codes 93922, 93923, and 93924 for billing and reimbursement purposes.

To obtain expert guidance on the appropriateness of these codes for specific systems, we highly recommend consulting with one of our experienced vascular specialists.

