

# simpleABI-400CL Automated Cuff-Link System

Newman Medical



Primary Care



Podiatry



Wound Care



Community Health



Rapid ABI-Q



Single-Level



TBI w/PPG

CLINIC PLACEMENT

EXAMS PERFORMED



American-Made Quality



Includes labor, parts, and accidental breakage



Automatic Cuff Selector



Push-Button Remote



simpleABI Roll Stand



ABI-Q Rapid Screen



DICOM Compatible

## simpleABI<sup>®</sup> with CUFF-LINK<sup>™</sup>

The simpleABI-400CL Automated Cuff-Link System offers accurate and efficient PAD testing for confident, early diagnosis. With versatile testing options, streamlined workflows, and simplified reimbursement processes, simpleABI enhances convenience and safety for both clinicians and patients. Backed by a 2-year warranty, it's the smart choice for proactive PAD management and improved patient outcomes.

### ABI-400CL Exams

- ABI-Q Rapid Screen
- Single-Level
- TBI with PPG

### ABI-400CL Includes:

- simpleABI software
- 4-port Cuff-Link with remote
- 8 MHz Doppler probe
- DPPG probe
- Cuffs
  - 2 x 10cm EZ-Clean Vasocuffs
  - 2 x 12cm EZ-Clean Vasocuffs
  - 2 x 1.9cm disposable digit cuffs
- Adjustable-height roll stand
- User manual / Quick reference guide
- Two-Year Warranty
  - Includes all parts, labor, and accidental breakage



## simpleABI Vascular Systems: Quality You Can Trust

### Made in the USA

- High-quality components, manufactured in Colorado.
- Durable design for busy clinical settings.
- 2-year warranty with complete coverage.

### Newman Medical: Your Partner in Vascular Care

- Leaders in ultrasound diagnostics since the 1970s.
- Committed to innovation and providing exceptional vascular products.
- Focused on service and long-term customer relationships.

### First-Class Support

- Experienced customer service team with extensive industry knowledge.
- Dedicated field representatives for personalized assistance.
- Comprehensive clinical support resources, both online and in-person.

MADE IN  
USA


## ABI-400CL



### Exams Performed



ABI-Q Rapid Screen



Single-Level



TBI with PPG

## Ordering Information

### System

#### SKU

ABI-400CL

#### DESCRIPTION

simpleABI Single-Level Cuff-Link System

### Accessories

#### SKU

ABI-120

STND-500

ACC-190

D5

SPPG

SFTW-100

ABI-WARR

#### DESCRIPTION

Microsoft® Windows® PC for use with simpleABI Systems

Variable-Height Roll Stand for Cuff-Link Systems

Large simpleABI carry bag

5Mhz vascular probe only

PPG Probe for use with Cuff-Link Systems - includes toe cuff

DICOM software integration

1 year additional warranty for simpleABI systems

### Replacement Cuffs

#### SKU

CUFF-200

CUFF-210

CUFF-210L

CUFF-220D

CUFF-240D

CUFF-SET-242

#### DESCRIPTION

EZ Clean Vascular Cuff, 10cm (ankle, small arm)

EZ Clean Vascular Cuff, 12cm (arm, ankle, calf)

EZ Clean Vascular Cuff, 12cm long (above knee)

Digit Cuff, 2.5cm - Latex-Free

Digit Cuff, 1.9cm - Latex-Free

EZ Clean Vascular Cuff Set - 2x10cm, 4x12cm, 2x12cm - 8/pk



## ABI-400CL Technical Specifications

**Height:** 107 cm (42 inches)

**Weight:** 5.5 kgs (12 pounds)

**Level of Protection against electrical shock:** Type B Applied Part, Class II Equipment

**Comply with the following Standards:** IEC60101-1, IEC60601-2, IEC60601-2-37

**Operating Temperature:** 10° ~ 40°C (50° ~ 104°F)

**Operating Humidity:** 30% ~ 75%

**Transport/Storage Temperature:** -20° ~ 50°C (-4° ~ 122°F)

**Transport/Storage Humidity:** 5% ~ 90%, non-condensing

**Battery Life:** 10 hours

**GTIN:** 00852785008630

## REIMBURSEMENT INFORMATION

Newman Medical's simpleABI vascular systems frequently utilize CPT codes 93922, 93923, and 93924 for billing and reimbursement purposes.

To obtain expert guidance on the appropriateness of these codes for specific systems, we highly recommend consulting with one of our experienced vascular specialists.

