

simpleABI-450CL Automated Cuff-Link System

Newman Medical



Primary Care



Podiatry



Wound Care



Community Health

CLINIC PLACEMENT



Rapid ABI-Q



Single-Level



TBI w/PPG



Exercise

EXAMS PERFORMED



American-Made Quality



Includes labor, parts, and accidental breakage



Automatic Cuff Selector



Push-Button Remote



simpleABI Roll Stand



ABI-Q Rapid Screen



DICOM Compatible

simpleABI[®] with CUFF-LINK[™]

The simpleABI-450CL Automated Cuff-Link System offers accurate and efficient PAD testing for confident, early diagnosis. With versatile testing options, streamlined workflows, and simplified reimbursement processes, simpleABI enhances convenience and safety for both clinicians and patients. Backed by a 2-year warranty, it's the smart choice for proactive PAD management and improved patient outcomes.

ABI-450CL Exams

- ABI-Q Rapid Screen
- Single-Level
- TBI with PPG
- Exercise

ABI-450CL Includes:

- simpleABI software
- 4-port Cuff-Link with remote
- 8 MHz Doppler probe
- DPPG probe
- Cuffs
 - 2 x 10cm EZ-Clean Vasocuffs
 - 2 x 12cm EZ-Clean Vasocuffs
 - 2 x 1.9cm disposable digit cuffs
- Adjustable-height roll stand
- User manual / Quick reference guide
- Two-Year Warranty
 - Includes all parts, labor, and accidental breakage



simpleABI Vascular Systems: Quality You Can Trust

Made in the USA

- High-quality components, manufactured in Colorado.
- Durable design for busy clinical settings.
- 2-year warranty with complete coverage.

Newman Medical: Your Partner in Vascular Care

- Leaders in ultrasound diagnostics since the 1970s.
- Committed to innovation and providing exceptional vascular products.
- Focused on service and long-term customer relationships.

First-Class Support

- Experienced customer service team with extensive industry knowledge.
- Dedicated field representatives for personalized assistance.
- Comprehensive clinical support resources, both online and in-person.




ABI-450CL



Exams Performed



ABI-Q Rapid Screen



Single-Level



TBI with PPG



Exercise Exams

Ordering Information

System

SKU

ABI-450CL

DESCRIPTION

simpleABI Single-Level Cuff-Link System with Exercise

Accessories

SKU

ABI-120

STND-500

ACC-190

D5

SPPG

SFTW-100

ABI-WARR

DESCRIPTION

Microsoft® Windows® PC for use with simpleABI Systems

Variable-Height Roll Stand for Cuff-Link Systems

Large simpleABI carry bag

5Mhz vascular probe only

PPG Probe for use with Cuff-Link Systems - includes toe cuff

DICOM software integration

1 year additional warranty for simpleABI systems

Replacement Cuffs

SKU

CUFF-200

CUFF-210

CUFF-210L

CUFF-220D

CUFF-240D

CUFF-SET-242

DESCRIPTION

EZ Clean Vascular Cuff, 10cm (ankle, small arm)

EZ Clean Vascular Cuff, 12cm (arm, ankle, calf)

EZ Clean Vascular Cuff, 12cm long (above knee)

Digit Cuff, 2.5cm - Latex-Free

Digit Cuff, 1.9cm - Latex-Free

EZ Clean Vascular Cuff Set - 2x10cm, 4x12cm, 2x12cm - 8/pk



ABI-450CL Technical Specifications

Height: 107 cm (42 inches)

Weight: 5.5 kgs (12 pounds)

Level of Protection against electrical shock: Type B Applied Part, Class II Equipment

Comply with the following Standards: IEC60101-1, IEC60601-2, IEC60601-2-37

Operating Temperature: 10° ~ 40°C (50° ~ 104°F)

Operating Humidity: 30% ~ 75%

Transport/Storage Temperature: -20° ~ 50°C (-4° ~ 122°F)

Transport/Storage Humidity: 5% ~ 90%, non-condensing

Battery Life: 10 hours

GTIN: 00852785008623

REIMBURSEMENT INFORMATION

Newman Medical's simpleABI vascular systems frequently utilize CPT codes 93922, 93923, and 93924 for billing and reimbursement purposes.

To obtain expert guidance on the appropriateness of these codes for specific systems, we highly recommend consulting with one of our experienced vascular specialists.

