

simpleABI-500CL Automated Cuff-Link System

Newman Medical



Specialists



Cardiology



Vascular Surgery



Interventional Radiology

CLINIC PLACEMENT



Rapid ABI-Q



Single-Level



Multi-Level



TBI w/PPG

EXAMS PERFORMED



American-Made Quality



Includes labor, parts, and accidental breakage



Automatic Cuff Selector



Push-Button Remote



simpleABI Roll Stand



ABI-Q Rapid Screen



DICOM Compatible

simpleABI with CUFF-LINK

The simpleABI-500CL Automated Cuff-Link System offers accurate and efficient PAD testing for confident, early diagnosis. With versatile testing options, streamlined workflows, and simplified reimbursement processes, simpleABI enhances convenience and safety for both clinicians and patients. Backed by a 2-year warranty, it's the smart choice for proactive PAD management and improved patient outcomes.

ABI-500CL Exams

- ABI-Q Rapid Screen
- Single-Level
- Multi-Level
- TBI with PPG

ABI-500CL Includes:

- simpleABI software
- 4-port Cuff-Link with remote
- 8 MHz Doppler probe
- DPPG probe
- Cuffs
 - 2 x 10cm EZ-Clean Vasocuffs
 - 4 x 12cm EZ-Clean Vasocuffs
 - 2 x 12cm long EZ-Clean Vasocuffs
 - 2 x 1.9cm disposable digit cuffs
- Adjustable-height roll stand
- User manual / Quick reference guide
- Two-Year Warranty
 - Includes all parts, labor, and accidental breakage



simpleABI Vascular Systems: Quality You Can Trust

Made in the USA

- High-quality components, manufactured in Colorado.
- Durable design for busy clinical settings.
- 2-year warranty with complete coverage.

Newman Medical: Your Partner in Vascular Care

- Leaders in ultrasound diagnostics since the 1970s.
- Committed to innovation and providing exceptional vascular products.
- Focused on service and long-term customer relationships.

First-Class Support

- Experienced customer service team with extensive industry knowledge.
- Dedicated field representatives for personalized assistance.
- Comprehensive clinical support resources, both online and in-person.

SIMPLIFY VASCULAR ASSESSMENTS WITH THE STEADYDOP PROBE

NEW Its unique stethoscope-like grip and optimized design ensure easy vessel tracking. Compatible with Newman Medical's DigiDop & simpleABI systems. Optional accessory. Contact your Newman Medical rep for details.




ABI-500CL



Exams Performed



ABI-Q Rapid Screen



Single-Level



Segmental



TBI with PPG

Ordering Information

System

SKU	DESCRIPTION
ABI-500CL	simpleABI Multi-Level Cuff-Link System

Cuff-Link Accessory Options

SKU	DESCRIPTION
ABI-120	Microsoft® Windows® PC for use with simpleABI Systems
STND-500	Variable-Height Roll Stand for Cuff-Link Systems
ACC-190	Large simpleABI carry bag
D5	5 Mhz vascular probe only
D5-SD	5 MHz SteadyDop vascular probe only
D8	8 MHz vascular probe only
D8-SD	8 MHz SteadyDop vascular probe only
SPPG	PPG Probe for use with Cuff-Link Systems - includes toe cuff
SFTW-100	DICOM software integration
ABI-WARR	1 year additional warranty for simpleABI systems

Replacement Cuffs

SKU	DESCRIPTION
CUFF-200	EZ Clean Vascular Cuff, 10cm (ankle, small arm)
CUFF-210	EZ Clean Vascular Cuff, 12cm (arm, ankle, calf)
CUFF-210L	EZ Clean Vascular Cuff, 12cm long (above knee)
CUFF-230	EZ-Clean Vascular Cuff, 17cm Contoured (thigh)
CUFF-250	EZ-Clean Vascular Cuff, 7cm long (metatarsal)
CUFF-220D	Digit Cuff, 2.5cm - Latex-Free
CUFF-240D	Digit Cuff, 1.9cm - Latex-Free
CUFF-SET-242	EZ Clean Vascular Cuff Set - 2x10cm, 4x12cm, 2x12cm - 8/pk



SIMPLIFY VASCULAR ASSESSMENTS WITH THE STEADYDOP PROBE

NEW Its unique stethoscope-like grip and optimized design ensure easy vessel tracking. Compatible with Newman Medical's DigiDop & simpleABI systems. Optional accessory. Contact your Newman Medical rep for details.



ABI-500CL Technical Specifications

Height:	107 cm (42 inches)
Weight:	5.5 kgs (12 pounds)
Level of Protection against electrical shock:	Type B Applied Part, Class II Equipment
Comply with the following Standards:	IEC60101-1, IEC60601-2, IEC60601-2-37
Operating Temperature:	10° ~ 40°C (50° ~ 104°F)
Operating Humidity:	30% ~ 75%
Transport/Storage Temperature:	-20° ~ 50°C (-4° ~ 122°F)
Transport/Storage Humidity:	5% ~ 90%, non-condensing
Battery Life:	10 hours
GTIN:	00852785008616

REIMBURSEMENT INFORMATION

Newman Medical's simpleABI vascular systems frequently utilize CPT codes 93922, 93923, and 93924 for billing and reimbursement purposes.

To obtain expert guidance on the appropriateness of these codes for specific systems, we highly recommend consulting with one of our experienced vascular specialists.

