



NO SURPRISES, JUST SUPPORT: DIGIDOP'S 5-YEAR WARRANTY

DigiDop 5-YEAR WARRANTY

Unexpected downtime can disrupt your workflow and impact patient care. That's why DigiDop offers an industry-leading 5-year warranty designed to get you back on track quickly and effortlessly.

The DigiDop 5-Year Warranty has you covered:

- **Complete Protection:** 5 years of coverage for all parts and labor—no extra costs.
- **Fast Turnaround:** Get back to your patients quickly with minimal disruption.
- **Zero Hassle:** No hidden fees, just easy and efficient service.
- **Dedicated Support:** We're here to help every step of the way.

Choose peace of mind. Choose DigiDop.



GO GREEN!

Save time and money while reducing your impact on the environment.

CHOOSE A RECHARGEABLE DIGIDOP



WHY CHOOSE DIGIDOP?

BUILT TO LAST

Durable, reliable, and backed by a comprehensive warranty. Made in the USA.

EXPERIENCE MATTERS

Decades of experience delivering innovative technology you can trust.

LEADING INNOVATION

Advanced features for crystal-clear sound and effortless use.

SMART & SUSTAINABLE

Rechargeable and long-lasting, saving you time and money.





Professional and Sensitive Dopplers



Over **45 years in the marketplace**; a testament to our longevity and success



Made in the USA for superior quality and performance



5-year all-inclusive warranty for maximum protection and minimized downtime, headaches



Best-in-class customer support; we always answer the phone

DigiDop Classic

- Standard size
- Superior battery endurance compared to competitors



DigiDop II

- Compact size
- Comfortable palm-fitting design
- Equivalent battery endurance to competitors



DigiDop 901

- Ideal for hospital use
- Clear stereo sound from dual top speakers
- Two probe cradles
- Standard drop-on charging pyramid
- Optional wall mount recharger



Vascular Probes



Model D5
5MHz Vascular Probe



Model D8
8MHz Vascular Probe



Model DPPG
Audio TBI PPG Probe

Fetal Probes



Model D2
2MHz Fetal Probe



Model D2W
2MHz Waterproof Fetal Probe




Model D3
3MHz Fetal Probe



Model D3W
3MHz Waterproof Fetal Probe

Interchangeable probes can be used on any DigiDop

Vascular




- Great sensitivity
- Ergonomic shape easy to hold, place in sterile glove finger for use in sterile fields
 - **D8 – 8 MHz** - best for shallow vessels like radial, ankle, foot, digits, penile
 - **D5 – 5 MHz** - wider beam, deeper penetration. Good for brachial, carotid, femoral, deep veins

Fetal



- Great sensitivity
- Smooth shape, comfortable for patient, easy to hold and clean
 - **D3 – 3 MHz** - wide angle, ideal for early fetal hearts
 - **D3W** - a unique 3 MHz waterproof probe
 - **D2 – 2 MHz** - deeper penetration
 - **D2W – 2 MHz** - waterproof

DPPG



- Unique audio PPG probe for use on digits. Great for toe-brachial studies

How to build your Doppler

DOPPLER CODE

PROBE CODE

Doppler

DD-300	Basic
DD-301	Basic + Recharge
DD-700	HR Display
DD-701	HR Display + Recharge
DD-901	Tabletop + HR Display + Recharge

DD-330	Basic
DD-330R	Basic + Recharge
DD-770	HR Display
DD-770R	HR Display + Recharge

Probe

D2	2MHz fetal probe
D2W	2MHz waterproof fetal probe
D3	3MHz fetal probe
D3W	3MHz waterproof fetal probe
D5	5MHz vascular probe
D8	8MHz vascular probe
DPPG	Audio PPG probe - includes toe cuff