



## DIGIDOP: WHERE VALUE MEETS CONFIDENCE BACKED BY AN INDUSTRY-LEADING 5-YEAR WARRANTY

### EXPERIENCE THE DIGIDOP DIFFERENCE:

- Reliable: 45+ years of proven Doppler expertise.
- Durable: American-made quality with a 5-year warranty.
- Versatile: Ideal for a wide range of applications.
- Efficient: Streamline workflows and enhance patient care.
- Supported: Exceptional clinical and technical support.

### VALUE BEYOND COMPARE:

- Competitive pricing with high-value features.
- Fast US delivery from our Colorado facility.
- Comprehensive 5-year warranty with no hidden fees or charges.

### KEY APPLICATIONS

- Emergency Medicine
- Labor & Delivery
- Vascular
- Cardiology
- Pediatrics/NICU



### RECOMMENDED DIGIDOP: DD-901 TABLETOP DOPPLER



### 5-YEAR WARRANTY

- All parts & labor
- No diagnostic fees
- No hidden fees
- No down time

Scan or click  
to explore  
DigiDop Dopplers.





# Professional and Sensitive Dopplers

Newman Medical



Over **45 years** in the marketplace; a testament to our longevity and success



**Made in the USA** for superior quality and performance



**5-year all-inclusive warranty** for maximum protection and minimized downtime, headaches



**Best-in-class customer support**; we always answer the phone

## DigiDop Classic

- Standard size
- Superior battery endurance compared to competitors



## DigiDop II

- Compact size
- Comfortable palm-fitting design
- Equivalent battery endurance to competitors



## DigiDop 901

- Ideal for hospital use
- Clear stereo sound from dual top speakers
- Two probe cradles
- Standard drop-on charging pyramid
- Optional wall mount recharger



### Vascular Probes



**Model D5**  
5MHz Vascular Probe



**Model D8**  
8MHz Vascular Probe



**Model DPPG**  
Audio TBI PPG Probe

### Fetal Probes



**Model D2**  
2MHz Fetal Probe



**Model D2W**  
2MHz Waterproof Fetal Probe




**Model D3**  
3MHz Fetal Probe



**Model D3W**  
3MHz Waterproof Fetal Probe

Interchangeable probes can be used on any DigiDop

### Vascular



- Great sensitivity
- Ergonomic shape easy to hold, place in sterile glove finger for use in sterile fields
  - **D8 – 8 MHz** - best for shallow vessels like radial, ankle, foot, digits, penile
  - **D5 – 5 MHz** - wider beam, deeper penetration. Good for brachial, carotid, femoral, deep veins

### Fetal



- Great sensitivity
- Smooth shape, comfortable for patient, easy to hold and clean
  - **D3 – 3 MHz** - wide angle, ideal for early fetal hearts
  - **D3W** - a unique 3 MHz waterproof probe
  - **D2 – 2 MHz** - deeper penetration
  - **D2W** - 2 MHz - waterproof

### DPPG



- Unique audio PPG probe for use on digits. Great for toe-brachial studies

### How to build your Doppler

DOPPLER CODE

PROBE CODE

#### Doppler

DD-300	Basic
DD-301	Basic + Recharge
DD-700	HR Display
DD-701	HR Display + Recharge
DD-901	Tabletop + HR Display + Recharge

DD-330	Basic
DD-330R	Basic + Recharge
DD-770	HR Display
DD-770R	HR Display + Recharge

#### Probe

D2	2MHz fetal probe
D2W	2MHz waterproof fetal probe
D3	3MHz fetal probe
D3W	3MHz waterproof fetal probe
D5	5MHz vascular probe
D8	8MHz vascular probe
DPPG	Audio PPG probe - includes toe cuff