



## EVERY EMS/FIRE DEPARTMENT NEEDS A DIGIDOP DOPPLER

In the fast-paced world of emergency response, every second counts. DigiDop Dopplers are a vital tool for EMS and Fire departments, providing crucial information quickly and reliably.

### How DigiDop Helps You:

- **Trauma:** Verify blood flow in fractures, crush injuries, and wounds.
- **Circulation:** Quickly find pulses, even in low-pressure situations.
- **Emergencies:** Portable and lightweight for any scene.
- **Pediatric & Neonatal:** Easily assess circulation in the smallest patients.
- **MCI (Mass Casualty Incidents):** Rapidly triage patients by assessing perfusion status.

### Benefits:

- **Easy to Use:** Minimal training required.
- **Enhanced Efficiency:** Quick decisions in critical situations.

### BEST DIGIDOP DOPPLER FOR EMS / FIRE DEPARTMENTS: OUR TOP PICK



DigiDop II

- COMPACT SIZE - FITS IN POCKET
- COMFORTABLE PALM-FITTING DESIGN
- EQUIVALENT BATTERY ENDURANCE TO COMPETITORS

VASCULAR



FETAL



DigiDop II

Smart Recharge  
System

MADE IN  
  
USA

5-YEAR  
WARRANTY

ALL PARTS & LABOR  
NO DIAGNOSTIC FEES  
NO HIDDEN FEES  
NO DOWN TIME

Scan or click  
to explore  
DigiDop Dopplers.



NM0150E



# Professional and Sensitive Dopplers

Newman Medical



Over **45 years** in the marketplace; a testament to our longevity and success



**Made in the USA** for superior quality and performance



**5-year all-inclusive warranty** for maximum protection and minimized downtime, headaches



**Best-in-class customer support**; we always answer the phone

## DigiDop Classic

- Standard size
- Superior battery endurance compared to competitors



## DigiDop II

- Compact size
- Comfortable palm-fitting design
- Equivalent battery endurance to competitors



## DigiDop 901

- Ideal for hospital use
- Clear stereo sound from dual top speakers
- Two probe cradles
- Standard drop-on charging pyramid
- Optional wall mount recharger



### Vascular Probes



**Model D5**  
5MHz Vascular Probe



**Model D8**  
8MHz Vascular Probe



**Model DPPG**  
Audio TBI PPG Probe

### Fetal Probes



**Model D2**  
2MHz Fetal Probe



**Model D2W**  
2MHz Waterproof Fetal Probe




**Model D3**  
3MHz Fetal Probe



**Model D3W**  
3MHz Waterproof Fetal Probe

Interchangeable probes can be used on any DigiDop

### Vascular




- Great sensitivity
- Ergonomic shape easy to hold, place in sterile glove finger for use in sterile fields
  - **D8 - 8 MHz** - best for shallow vessels like radial, ankle, foot, digits, penile
  - **D5 - 5 MHz** - wider beam, deeper penetration. Good for brachial, carotid, femoral, deep veins

### Fetal



- Great sensitivity
- Smooth shape, comfortable for patient, easy to hold and clean
  - **D3 - 3 MHz** - wide angle, ideal for early fetal hearts
  - **D3W** - a unique 3 MHz waterproof probe
  - **D2 - 2 MHz** - deeper penetration
  - **D2W** - 2 MHz - waterproof

### DPPG



- Unique audio PPG probe for use on digits. Great for toe-brachial studies

### How to build your Doppler

DOPPLER CODE

PROBE CODE

#### Doppler

DD-300	Basic
DD-301	Basic + Recharge
DD-700	HR Display
DD-701	HR Display + Recharge
DD-901	Tabletop + HR Display + Recharge

DigiDop "Classic"

DigiDop II

DD-330	Basic
DD-330R	Basic + Recharge
DD-770	HR Display
DD-770R	HR Display + Recharge

#### Probe

D2	2MHz fetal probe
D2W	2MHz waterproof fetal probe
D3	3MHz fetal probe
D3W	3MHz waterproof fetal probe
D5	5MHz vascular probe
D8	8MHz vascular probe
DPPG	Audio PPG probe - includes toe cuff