


# natus ↔ DigiDop Cross-Reference

Learn how Newman Medical started and developed the Natus Doppler line.



	natus Model	REPLACEMENT Newman Medical DigiDop Classic	REPLACEMENT Newman Medical DigiDop II
NON-DISPLAY	<p><b>Elite 100</b></p>	<p>DD-300</p>	<p>DD-330</p>
	<p><b>Elite 100/100R</b></p>	<p>DD-301</p>	<p>DD-330R</p>
	<p><b>PD-II</b></p>	<p>DD-301</p>	<p>DD-330R</p>
DISPLAY	<p><b>Elite 200</b></p>	<p>DD-700</p>	<p>DD-770</p>
	<p><b>Elite 200/200R</b></p>	<p>DD-701</p>	<p>DD-770R</p>
TABLETOP	<p><b>IMEXDOP CT+</b></p>	<p>REPLACEMENT Newman Medical Model</p> <p>DD-901</p>	



Includes labor, parts, and accidental breakage



American-Made Quality



Includes cords and probes



Waterproof Option



Rechargeable Option

H A N D H E L D D O P L E R S N O N - D I S P L A Y

T A B L E T O P D I S P L A Y

P R O B E O R D E R I N G I N F O R M A T I O N O N B A C K S I D E

## Vascular Probes



**Model D5**  
5MHz Vascular Probe



**Model D8**  
8MHz Vascular Probe



**Model DPPG**  
Audio TBI PPG Probe

## Fetal Probes



**Model D2**  
2MHz Fetal Probe



**Model D2W**  
2MHz Waterproof Fetal Probe



**Model D3**  
3MHz Fetal Probe



**Model D3W**  
3MHz Waterproof Fetal Probe

## Professional and Sensitive Dopplers

### Quality

- Better sound - Digital Sound Clarity™ produces up to 6x signal of other Dopplers
- Better warranty - 5 year honest warranty all parts & labor. Includes rechargeable batteries on rechargeable units.
- Better support & service - knowledgeable staff and guaranteed service

### Sustainability

- Smart Recharge System™ (SRS)
  - Significant cost savings
  - New rechargeable batteries hold charge for years
  - May be recharged more than 2000 times

### Innovation

- Continual product development and research

### DigiDop Classic



- Standard size
- Longer battery life than similar Natus products

### DigiDop II



- Smaller size
- Ergonomic - fits comfortably in your hand
- Similar battery life to Natus

### DigiDop 901



- Ideal for hospital use
- Great room-filling sound from two speakers on top
- Two probe cradles
- Standard drop-on charging pyramid
- Optional wall mount recharger

## How to build your Doppler

DOPPLER CODE

PROBE CODE

### Doppler

DD-300	Basic
DD-301	Basic + Recharge
DD-700	HR Display
DD-701	HR Display + Recharge
DD-901	Tabletop + HR Display + Recharge

### Probe

D2	2MHz fetal probe
D2W	2MHz waterproof fetal probe
D3	3MHz fetal probe
D3W	3MHz waterproof fetal probe
D5	5MHz vascular probe
D8	8MHz vascular probe
DPPG	Audio PPG probe - includes toe cuff

## Interchangeable probes can be used on any DigiDop

### Vascular



- Great sensitivity
- Ergonomic shape easy to hold, place in sterile glove finger for use in sterile fields
- D8 - 8 MHz** - best for shallow vessels like radial, ankle, foot, digits, penile
- D5 - 5 MHz** - wider beam, deeper penetration. Good for brachial, carotid, femoral, deep veins

### Fetal



- Great sensitivity
- Smooth shape, comfortable for patient, easy to hold and clean
- D3 - 3 MHz** - wide angle, ideal for early fetal hearts
- D3W** - a unique 3 MHz waterproof probe
- D2 - 2MHz** - deeper penetration
- D2W** - 2 MHz - waterproof

### PPG



- DPPG**
- Unique audio PPG probe for use on digits. Great for toe-brachial studies

## NATUS-NEWMAN MEDICAL PART NUMBER CROSS-REFERENCE



**Elite**  
100/100R

Natus	DigiDop Classic	DigiDop II
EN20	DD-300-D2	DD-330-D2
EN20R	DD-301-D2	DD-330R-D2
EN30	DD-300-D3	DD-330-D3
EN30R	DD-301-D3	DD-330R-D3
EN50	DD-300-D5	DD-330-D5
EN50R	DD-301-D5	DD-330R-D5
EN58	DD-300-VASC	--
EN80	DD-300-D8	DD-330-D8
EN80R	DD-301-D8	DD-330R-D8
EN8A	DD-PAD	--



**Elite**  
200/200R

Natus	DigiDop Classic	DigiDop II
ED20	DD-700-D2	DD-770-D2
ED20R	DD-701-D2	DD-770R-D2
ED23	DD-700-OB	--
ED2W	DD-700-D2W	DD-770-D2W
ED2WR	DD-701-D2W	DD-770R-D2W
ED30	DD-700-D3	DD-770-D3
ED30R	DD-701-D3	DD-770R-D3
ED50	DD-700-D5	DD-770-D5
ED50R	DD-701-D5	DD-770R-D5
ED58	DD-700-VASC	--
ED80	DD-700-D8	DD-770-D8
ED80R	DD-701-D8	DD-770R-D8



**Pocket-Dop II**

Natus	DigiDop Classic	DigiDop II
P200	DD-301-D2	DD-330R-D2
P300	DD-301-D3	DD-330R-D3
P500	DD-301-D5	DD-330R-D5
P800	DD-301-D8	DD-330R-D8



**IMEXDOP CT+**

Natus	DigiDop Classic	DigiDop II
CT20+	DD-901-D2	--
CT30+	DD-901-D3	--
CT50+	DD-901-D5	--
CT80+	DD-901-D8	--

