

Pro-Lab Vascular Bundle



Specialists



Cardiology



Vascular Surgery



Interventional Radiology

CLINIC PLACEMENT



Rapid ABI-Q



Single-Level



Multi-Level



TBI w/PPG



Exercise

EXAMS PERFORMED



American-Made Quality



Includes labor, parts, and accidental breakage



Automatic Cuff Selector



Push-Button Remote



simpleABI Roll Stand



ABI-Q Rapid Screen



DICOM

simpleABI[®] with CUFF-LINK[™]

The Pro-Lab Bundle delivers the complete simpleABI-600CL testing platform – every exam type, every billing code – with one defining upgrade: a 5-year warranty in place of the standard two. For high-volume vascular labs, health systems, and practices where this system becomes mission-critical infrastructure, the extended coverage means uninterrupted clinical operation and long-term cost certainty. It's the same comprehensive PAD testing capability your team depends on, protected for the long haul.

Pro-Lab Exams

- ABI-Q Rapid Screen
- Single-Level
- Multi-Level
- TBI with PPG
- Exercise & Stress

Pro-Lab Includes:

- simpleABI software
- 17" touchscreen Microsoft® Windows® PRO PC
- 8-port Cuff-Link with remote
- 8 MHz Doppler pencil probe
- 5 MHz Doppler pencil probe
- 5 MHz SteadyDop ABI probe
- DPPG probe
- Cuffs
 - 2 x 10cm EZ-Clean Vasocuffs
 - 4 x 12cm EZ-Clean Vasocuffs
 - 2 x 12cm long EZ-Clean Vasocuffs
 - 2 x 17cm contoured EZ-Clean Vasocuffs
 - 2 x 7cm EZ-Clean Vasocuffs
 - 2 x 1.9cm disposable digit cuffs
 - 2 x 2.5cm disposable digit cuffs
- Adjustable-height roll stand
- Five-Year Warranty
 - Includes all parts, labor, and accidental breakage
- DICOM Software Integration
 - Full worklist capability
 - Image and structured reporting
- User manual / Quick reference guide



simpleABI Vascular Systems: Quality You Can Trust

Made in the USA

- High-quality components, manufactured in Colorado.
- Durable design for busy clinical settings.
- 2-year warranty with complete coverage.

Newman Medical: Your Partner in Vascular Care

- Leaders in ultrasound diagnostics since the 1970s.
- Committed to innovation and providing exceptional vascular products.
- Focused on service and long-term customer relationships.

First-Class Support

- Experienced customer service team with extensive industry knowledge.
- Dedicated field representatives for personalized assistance.
- Comprehensive clinical support resources, both online and in-person.




ABI-600CL




Exams Performed



ABI-Q Rapid Screen



Single-Level



Segmental



TBI with PPG



Exercise Exams

Ordering Information

System

SKU
PRO-LAB

DESCRIPTION

Pro-Lab Multi-Level with Exercise Cuff-Link System

PRO-LAB Accessory Options

SKU

ABI-120
STND-500
ACC-190
D5
D5-SD
D8
D8-SD
SPPG
SFTW-100

DESCRIPTION

Microsoft® Windows® PC for use with simpleABI Systems
Variable-Height Roll Stand for Cuff-Link Systems
Large simpleABI carry bag
5 Mhz vascular probe only
5 MHz SteadyDop vascular probe only
8 MHz vascular probe only
8 MHz SteadyDop vascular probe only
PPG Probe for use with Cuff-Link Systems - includes toe cuff
DICOM software integration



Replacement Cuffs

SKU

CUFF-200
CUFF-210
CUFF-210L
CUFF-230
CUFF-250
CUFF-220D
CUFF-240D
CUFF-SET-242

DESCRIPTION

EZ Clean Vascular Cuff, 10cm (ankle, small arm)
EZ Clean Vascular Cuff, 12cm (arm, ankle, calf)
EZ Clean Vascular Cuff, 12cm long (above knee)
EZ-Clean Vascular Cuff, 17cm Contoured (thigh)
EZ-Clean Vascular Cuff, 7cm long (metatarsal)
Digit Cuff, 2.5cm - Latex-Free
Digit Cuff, 1.9cm - Latex-Free
EZ Clean Vascular Cuff Set - 2x10cm, 4x12cm, 2x12cm - 8/pk



SIMPLIFY VASCULAR ASSESSMENTS WITH THE STEADYDOP PROBE

Its unique stethoscope-like grip and optimized design ensure easy vessel tracking. Compatible with Newman Medical's DigiDop & simpleABI systems. Optional accessory. Contact your Newman Medical rep for details.



ABI-600CL Technical Specifications

Height: 107 cm (42 inches)

Weight: 5.5 kgs (12 pounds)

Level of Protection against electrical shock: Type B Applied Part, Class II Equipment

Comply with the following Standards: IEC60101-1, IEC60601-2, IEC60601-2-37

Operating Temperature: 10° ~ 40°C (50° ~ 104°F)

Operating Humidity: 30% ~ 75%

Transport/Storage Temperature: -20° ~ 50°C (-4° ~ 122°F)

Transport/Storage Humidity: 5% ~ 90%, non-condensing

Battery Life: 10 hours

GTIN: 852785008609

REIMBURSEMENT INFORMATION

Newman Medical's simpleABI vascular systems frequently utilize CPT codes 93922, 93923, and 93924 for billing and reimbursement purposes.

To obtain expert guidance on the appropriateness of these codes for specific systems, we highly recommend consulting with one of our experienced vascular specialists.

