

# simpleABI<sup>300</sup>

## Manual Doppler-Based Single-Level ABI + TBI Testing

Newman Medical

**2 YEAR WARRANTY**  
ALL PARTS & LABOR  
ACCIDENTAL BREAKAGE  
NO DIAGNOSTIC FEES  
NO HIDDEN FEES  
MINIMAL DOWN TIME

MADE IN  
  
USA

**REAL PEOPLE  
REAL SUPPORT**

NO ROBOTS.  
NO RUNAROUND.  
JUST REAL HELP.



The simpleABI-300 manual system is an excellent choice for budget-conscious practices that require more advanced PAD testing capabilities. This versatile system allows clinicians to perform both single-level ABI exams and Toe Brachial Index (TBI) measurements, providing a comprehensive assessment of lower extremity arterial health. Its affordability and portability make it an ideal solution for a variety of clinical settings.

It's backed by the simpleABI 2-year warranty that **covers everything** from parts to labor, with no hidden fees or diagnostic charges.


### Expand Your Diagnostic Abilities Without Breaking the Bank

#### simpleABI-300 features:

- Made in the USA with high-quality materials
- Comprehensive 2-year warranty coverage
- Great support: on-site, online or by phone
- Comprehensive PAD testing with ABI and TBI.
- Accurate TBI results, even with calcified arteries.
- Small vessel assessment via toe blood flow.
- Affordable, advanced PAD testing.
- Portable and compact design.



#### simpleABI-300 Exams



Doppler ABI



TBI

- Ideal for:**
- Primary care clinics
  - Podiatry offices
  - Wound care centers
  - Vascular clinics
  - Budget-conscious practices needing comprehensive PAD testing (ABI & TBI)

#### simpleABI-300 Includes:

- simpleABI software
- Doppler with 8 MHz probe
- DPPG probe
- PVR
- Manual cuff inflator
- Cuffs
  - 2 x 10cm EZ-Clean Vasocuffs
  - 2 x 12cm EZ-Clean Vasocuffs
  - 1 x 1.9cm digit cuff
- Carry bag
- User manual / Quick reference guide
- Two-year warranty
  - Includes all parts, labor, and accidental breakage

#### simpleABI-300 Accessories:

SKU	Description
ABI-120	Microsoft® Windows® PC
D5	5 MHz vascular probe only
D8	8 MHz vascular probe only
STND-130	simpleABI roll stand for manual systems
ACC-175	simpleABI-300 carry bag
PUMP-100	Manual cuff inflator

#### simpleABI-300 order info:

SKU	Description
ABI-300	simpleABI-300 manual system