

simpleABI[®] 400CL

with CUFF-LINK™

Automated Doppler-Based Single-Level ABI + TBI Testing

Dual-Testing Capabilities: Includes Doppler-Free PVR-Based Rapid-Test

INDUSTRY-LEADING YEAR WARRANTY

2

ALL PARTS & LABOR
ACCIDENTAL BREAKAGE
NO DIAGNOSTIC FEES
NO HIDDEN FEES
MINIMAL DOWN TIME

MADE IN
USA

REAL PEOPLE REAL SUPPORT

NO ROBOTS.
NO RUNAROUND.
JUST REAL HELP.



The simpleABI-400CL automated Cuff-Link system is a sophisticated solution for clinics seeking efficient PAD testing. Offering dual testing capabilities, it performs both the Doppler-free, rapid ABI-Q exam (allowing any staff member to quickly screen patients) and comprehensive, Doppler-based ABI and TBI exams. Integrating automated cuff selection and inflation with these advanced capabilities streamlines workflows, minimizes manual steps, and boosts diagnostic accuracy, freeing clinicians to focus on patient care and precise interpretation.

simpleABI-400CL Features

- Made in the USA with high-quality materials
- Comprehensive 2-year warranty coverage
- Expert support: on-site, online and phone assistance
- Dual testing: ABI-Q rapid-test and Doppler-based exams.
- Automated cuff selection and inflation to reduce user error.
- Versatile testing capabilities for various patient needs.
- Accurate and reliable results.
- User-friendly interface with minimal training required.
- Streamlined workflow

Ideal for:

- Primary care clinics
- Podiatry offices
- Wound care centers
- Community Health



simpleABI-400CL Exams



Rapid ABI-Q exam



ABI-Q plus exercise
(no gym equipment needed)



Doppler ABI



TBI

simpleABI-400CL Includes:

- simpleABI software
- 4-port Cuff-Link with remote
- 8 MHz Doppler probe
- DPPG probe
- Cuffs
 - 2 x 10cm EZ-Clean Vasocuffs
 - 2 x 12cm EZ-Clean Vasocuffs
 - 2 x 1.9cm disposable digit cuffs
- Adjustable-height roll stand
- User manual / Quick reference guide
- Two-year warranty
 - Includes all parts, labor, and accidental breakage

simpleABI-400CL Accessories:

SKU

ABI-120
STND-500
ACC-190
D5
D8
SPPG
SFTW-100
ABI-WARR

Description

Microsoft® Windows® PC
Variable-height roll stand for Cuff-Link systems
Large simpleABI carry bag
5MHz vascular probe only
8 MHz vascular probe only
Cuff-Link PPG probe - includes toe cuff
DICOM software integration
1 year additional warranty for simpleABI systems

simpleABI-400CL Order Info:

SKU	Description
ABI-400CL	simpleABI-400CL Cuff-Link system