

simpleABI[®] 400CL

with CUFF-LINK[™]

Automated Doppler-Based Single-Level ABI + TBI Testing

Dual-Testing Capabilities: Includes Doppler-Free PVR-Based Rapid-Test

Newman Medical

INDUSTRY-LEADING 2 YEAR WARRANTY

ALL PARTS & LABOR
ACCIDENTAL BREAKAGE
NO DIAGNOSTIC FEES
NO HIDDEN FEES
MINIMAL DOWN TIME

MADE IN
USA

REAL PEOPLE REAL SUPPORT

NO ROBOTS.
NO RUNAROUND.
JUST REAL HELP.



The simpleABI-400CL automated Cuff-Link system is a sophisticated solution for clinics seeking efficient PAD testing. Offering dual testing capabilities, it performs both the Doppler-free, rapid ABI-Q exam (allowing any staff member to quickly screen patients) and comprehensive, Doppler-based ABI and TBI exams. Integrating automated cuff selection and inflation with these advanced capabilities streamlines workflows, minimizes manual steps, and boosts diagnostic accuracy, freeing clinicians to focus on patient care and precise interpretation.

simpleABI-400CL Features



- Made in the USA with high-quality materials
- Comprehensive 2-year warranty coverage
- Expert support: on-site, online and phone assistance
- Dual testing: ABI-Q rapid-test and Doppler-based exams.
- Automated cuff selection and inflation to reduce user error.
- Versatile testing capabilities for various patient needs.
- Accurate and reliable results.
- User-friendly interface with minimal training required.
- Streamlined workflow

Ideal for:

- Primary care clinics
- Podiatry offices
- Wound care centers
- Community Health



simpleABI-400CL Exams



simpleABI-400CL Includes:

- simpleABI software
- 4-port Cuff-Link with remote
- 8 MHz Doppler probe
- SPPG probe
- Cuffs
 - 2 x 10cm EZ-Clean Vasocuffs
 - 2 x 12cm EZ-Clean Vasocuffs
 - 2 x 1.9cm disposable digit cuffs
- Adjustable-height roll stand
- User manual / Quick reference guide
- Two-year warranty
 - Includes all parts, labor, and accidental breakage

simpleABI-400CL Accessories:

SKU	Description
ABI-120	Microsoft® Windows® PC
STND-500	Variable-height roll stand for Cuff-Link systems
ACC-190	Large simpleABI carry bag
D5	5 MHz vascular probe only
D5-SD	5 MHz SteadyDop vascular probe only
D8	8 MHz vascular probe only
D8-SD	8 MHz SteadyDop vascular probe only
SPPG	Cuff-Link PPG probe - includes toe cuff
SFTW-100	DICOM software integration
ABI-WARR	1 year additional warranty for simpleABI systems

simpleABI-400CL Order Info:

SKU	Description
ABI-400CL	simpleABI-400CL Cuff-Link system



SIMPLIFY VASCULAR ASSESSMENTS WITH THE STEADYDOP PROBE

Its unique stethoscope-like grip and optimized design ensure easy vessel tracking. Compatible with Newman Medical's DigiDop & simpleABI systems. Optional accessory. Contact your Newman Medical rep for details.

