

# simpleABI<sup>with CUFF-LINK</sup> 450CL

## Automated Doppler-Based Single-Level ABI + TBI Testing with Stress/Exercise

**Dual-Testing Capabilities:** Includes Doppler-Free PVR-Based Rapid-Test

### INDUSTRY-LEADING YEAR WARRANTY

2

ALL PARTS & LABOR  
ACCIDENTAL BREAKAGE  
NO DIAGNOSTIC FEES  
NO HIDDEN FEES  
MINIMAL DOWN TIME



### REAL PEOPLE REAL SUPPORT

NO ROBOTS.  
NO RUNAROUND.  
JUST REAL HELP.



The simpleABI-450CL automated Cuff-Link system offers a versatile and efficient solution for comprehensive PAD assessment, going beyond basic ABI testing. Combining the convenience of automated cuff inflation with a broad range of exam capabilities, including exercise and stress testing, this system excels at uncovering hidden PAD and providing more accurate diagnoses. Its dual testing capability includes the rapid, Doppler-free ABI-Q for quick screening by any staff member, as well as comprehensive Doppler-based exams for in-depth analysis.

### simpleABI-450CL Features

- Made in the USA with high-quality materials
- Comprehensive 2-year warranty coverage
- Expert support: on-site, online and phone assistance
- Dual testing: ABI-Q rapid-test and Doppler-based exams.
- Comprehensive PAD testing with exercise/stress testing.
- Automated cuff selection and inflation to reduce user error.
- Accurate diagnoses.
- Improved workflow.

### Ideal for:

- Primary care clinics
- Podiatry offices
- Wound care centers
- Community Health




### simpleABI-450CL Exams



Rapid ABI-Q exam



ABI-Q plus exercise  
(no gym equipment needed)



Doppler ABI



TBI



Stress/Exercise

### simpleABI-450CL Includes:

- simpleABI software
- 4-port Cuff-Link with remote
- 8 MHz Doppler probe
- SPPG probe
- Cuffs
  - 2 x 10cm EZ-Clean Vasocuffs
  - 2 x 12cm EZ-Clean Vasocuffs
  - 2 x 1.9cm disposable digit cuffs
- Adjustable-height roll stand
- User manual / Quick reference guide
- Two-year warranty
  - Includes all parts, labor, and accidental breakage

### simpleABI-450CL Accessories:

#### SKU

ABI-120  
STND-500  
ACC-190  
D5  
D8  
SPPG  
SFTW-100  
ABI-WARR

#### Description

Microsoft® Windows® PC  
Variable-height roll stand for Cuff-Link systems  
Large simpleABI carry bag  
5MHz vascular probe only  
8 MHz vascular probe only  
Cuff-Link PPG probe - includes toe cuff  
DICOM software integration  
1 year additional warranty for simpleABI systems

### simpleABI-450CL Order Info:

SKU	Description
ABI-450CL	simpleABI-450CL Cuff-Link system