

simpleABI^{500CL} with CUFF-LINK™

Newman Medical

Automated Doppler-Based Segmental ABI + TBI Testing

Dual-Testing Capabilities: Includes Doppler-Free PVR-Based Rapid-Test

INDUSTRY-LEADING YEAR WARRANTY

2

ALL PARTS & LABOR
ACCIDENTAL BREAKAGE
NO DIAGNOSTIC FEES
NO HIDDEN FEES
MINIMAL DOWN TIME

MADE IN

USA

REAL PEOPLE REAL SUPPORT

NO ROBOTS.
NO RUNAROUND.
JUST REAL HELP.



The simpleABI-500CL automated Cuff-Link system is an advanced PAD testing solution designed for clinicians requiring precise arterial blockage localization. Combining efficient automated cuff inflation with comprehensive segmental exam capabilities, it provides detailed leg blood flow mapping for accurate diagnosis and targeted treatment planning. 1 Like all Cuff-Link systems, it also offers dual testing capabilities, including the rapid, Doppler-free ABI-Q for quick screening by any staff member, and in-depth Doppler-based exams.

simpleABI-500CL Features

- Made in the USA with high-quality materials
- Comprehensive 2-year warranty coverage
- Great support: on-site, online or by phone
- Dual testing: ABI-Q and Doppler exams.
- Advanced diagnostics with segmental exams to pinpoint blockages.
- Automated cuff selection and inflation at multiple levels to reduce user error.

Ideal for:

- Vascular Specialists
- Cardiology
- Vascular Surgery
- Interventional Radiology




simpleABI-500CL Exams




Rapid ABI-Q exam



ABI-Q plus exercise
(no gym equipment needed)



Doppler ABI



Segmental Pressures



TBI

simpleABI-500CL Includes:

- simpleABI software
- 8-port Cuff-Link with remote
- 8 MHz Doppler probe
- SPPG probe
- Cuffs
 - 2 x 10cm EZ-Clean Vasocuffs
 - 4 x 12cm EZ-Clean Vasocuffs
 - 2 x 12cm long EZ-Clean Vasocuffs
 - 2 x 1.9cm disposable digit cuffs
- Adjustable-height roll stand
- User manual / Quick reference guide
- Two-Year Warranty
 - Includes all parts, labor, and accidental breakage

simpleABI-500CL Accessories:

SKU

ABI-120
STND-500
ACC-190
D5
D8
SPPG
SFTW-100
ABI-WARR

Description

Microsoft® Windows® PC
Variable-height roll stand for Cuff-Link systems
Large simpleABI carry bag
5Mhz vascular probe only
8 MHz vascular probe only
Cuff-Link PPG probe - includes toe cuff
DICOM software integration
1 year additional warranty for simpleABI systems

simpleABI-500CL Order Info:

SKU	Description
ABI-500CL	simpleABI-500CL Cuff-Link system