

simpleABI[®] 600CL

with CUFF-LINK™

Newman Medical

Automated Doppler-Based Segmental ABI + TBI Testing with Stress/Exercise
Dual-Testing Capabilities: Includes Doppler-Free PVR-Based Rapid-Test

INDUSTRY-LEADING YEAR WARRANTY

2

ALL PARTS & LABOR
 ACCIDENTAL BREAKAGE
 NO DIAGNOSTIC FEES
 NO HIDDEN FEES
 MINIMAL DOWN TIME

MADE IN

USA

REAL PEOPLE REAL SUPPORT

NO ROBOTS.
 NO RUNAROUND.
 JUST REAL HELP.



The simpleABI-600CL Automated Cuff-Link system represents the pinnacle of PAD testing technology. This comprehensive system combines automated efficiency with the widest range of simpleABI exams, including ABI-Q, single-level, segmental, TBI, stress, and exercise testing. This breadth of options empowers clinicians to conduct thorough PAD evaluations, uncover hidden disease, pinpoint blockage locations, and optimize treatment plans with unparalleled accuracy and efficiency. And like all Cuff-Link systems, it offers dual testing capabilities, including the rapid, Doppler-free ABI-Q for quick screening by any staff member, and in-depth Doppler-based exams.

simpleABI-600CL Features


- Made in the USA with high-quality materials
- Comprehensive 2-year warranty coverage
- Great support: on-site, online or by phone
- Dual testing: ABI-Q and Doppler exams.
- All PAD tests, including ABI-Q, Doppler ABI with exercise/stress test, TBI, and segmental exams.
- Automated cuff selection and inflation at all levels to reduce user error.

Ideal for:

- Vascular specialists
- Cardiology
- Vascular surgery
- Interventional radiology



simpleABI-600CL Exams



simpleABI-600CL Includes:

- simpleABI software
- 8-port Cuff-Link with remote
- 8 MHz Doppler probe
- SPPG probe
- Cuffs
 - 2 x 10cm EZ-Clean Vasocuffs
 - 4 x 12cm EZ-Clean Vasocuffs
 - 2 x 12cm long EZ-Clean Vasocuffs
 - 2 x 1.9cm disposable digit cuffs
- Adjustable-height roll stand
- User manual / Quick reference guide
- Two-year warranty
 - Includes all parts, labor, and accidental breakage

simpleABI-600CL Accessories:

| SKU | Description |
|----------|--|
| ABI-120 | Microsoft® Windows® PC |
| STND-500 | Variable-height roll stand for Cuff-Link systems |
| ACC-190 | Large simpleABI carry bag |
| D5 | 5 Mhz vascular probe only |
| D8 | 8 MHz vascular probe only |
| SPPG | Cuff-Link PPG probe - includes toe cuff |
| SFTW-100 | DICOM software integration |
| ABI-WARR | 1 year additional warranty for simpleABI systems |

simpleABI-600CL Order Info:

| SKU | Description |
|-----------|----------------------------------|
| ABI-600CL | simpleABI-600CL Cuff-Link system |