

simpleABI^{with CUFF-LINK™} PRO-LAB BUNDLE

Newman Medical

Automated Doppler-Based Segmental ABI + TBI Testing with Stress/Exercise
Dual-Testing Capabilities: Includes Doppler-Free PVR-Based Rapid-Test

INDUSTRY-LEADING 5 YEAR WARRANTY

ALL PARTS & LABOR
 ACCIDENTAL BREAKAGE
 NO DIAGNOSTIC FEES
 NO HIDDEN FEES
 MINIMAL DOWN TIME



REAL PEOPLE REAL SUPPORT

NO ROBOTS.
 NO RUNAROUND.
 JUST REAL HELP.



The Pro-Lab Bundle delivers the ultimate simpleABI-600CL testing platform—covering every exam type and billing code—with one defining upgrade: a 5-year warranty that includes labor, parts, and accidental breakage. Designed for high-volume vascular labs and health systems where equipment is mission-critical, this bundle ensures uninterrupted clinical operation and long-term cost certainty. It provides the same hospital-grade performance and comprehensive PAD testing capability your team depends on, protected for the long haul at less than 50% of the cost of traditional vascular lab systems.

PRO-LAB Bundle Features

- Made in the USA with high-quality materials
- Comprehensive 5-year warranty coverage
- Great support: on-site, online or by phone
- Dual testing: ABI-Q and Doppler exams.
- All PAD tests, including ABI-Q, Doppler ABI with exercise/stress test, TBI, and segmental exams.
- Automated cuff selection and inflation at all levels to reduce user error.
- DICOM software integration with full worklist capability.

Ideal for:

- Vascular labs
- Vascular specialists
- Cardiology
- Vascular surgery
- Interventional radiology



Integrated DICOM
 5-Year Warranty
 Windows PRO PC

Pro-Lab Vascular Bundle Exams



Rapid ABI-Q exam



ABI-Q plus exercise
 (no gym equipment needed)



Doppler ABI



Segmental Pressures



TBI



Stress/Exercise

Pro-Lab Bundle Includes:

- simpleABI software
- 17" touchscreen Microsoft® Windows® PRO PC
- 8-port Cuff-Link with remote
- 8 MHz Doppler pencil probe
- 5 MHz Doppler pencil probe
- 5 MHz SteadyDop ABI probe
- DPPG probe
- Cuffs
 - 2 x 10cm EZ-Clean Vasocuffs
 - 4 x 12cm EZ-Clean Vasocuffs
 - 2 x 12cm long EZ-Clean Vasocuffs
 - 2 x 17cm contoured EZ-Clean Vasocuffs
 - 2 x 7cm EZ-Clean Vasocuffs
 - 2 x 1.9cm disposable digit cuffs
 - 2 x 2.5cm disposable digit cuffs
- Adjustable-height roll stand
- Five-Year Warranty
 - Includes all parts, labor, and accidental breakage
- DICOM Software Integration
 - Full worklist capability
 - Image and structured reporting
- User manual / Quick reference guide

Pro-Lab Bundle Accessories:

SKU	Description
ABI-120	Microsoft® Windows® PC
STND-500	Variable-height roll stand for Cuff-Link systems
ACC-190	Large simpleABI carry bag
D5	5 MHz vascular probe only
D5-SD	5 MHz SteadyDop vascular probe only
D8	8 MHz vascular probe only
D8-SD	8 MHz SteadyDop vascular probe only
SPPG	Cuff-Link PPG probe - includes toe cuff
SFTW-100	DICOM software integration
ABI-WARR	1 year additional warranty for simpleABI systems

Pro-Lab Bundle Order Info:

SKU	Description
PRO-LAB	Pro-Lab Complete Vascular Bundle