

File Menu

Beginning a new exam

From the computer Desktop, select the simpleABI icon. The simpleABI start-up screen with menu items on the top will appear. Select File→New and the exam to be performed. Depending on the system chosen, you will have different choices of exams. Choose an exam and it will open on the screen.

Loading an Exam

From the simpleABI start-up screen, Select File→Load and the exam to be performed. You may select any previous report to populate the loading exam with data relevant to the exam. I.e., you can load an ABI exam into a segmental exam; all patient information and arm/ankle pressures and waveforms from previous exam will be loaded. You would need only to perform upper leg pressure/waveform measurements. Depending on the system chosen, you will have different choices of exams. Choose an exam to load, select a previous exam and it will open on the screen.

Saving an Exam

You can save your finished report using CTRL+S or selecting File→Save. Files are appended with Patient ID and a timestamp to differentiate exams. You may also overwrite this default file name. By default, the save command automatically saves the exam in PDF and XML data files. You may select additional file types under Edit→Output Formats (TIFF, BMP, JPG, CSV)

Printing an Exam

You can print your exams using CTRL+P or selecting File→Print

Closing an exam

You can close any open exam or editing form by selecting File→Close Form or CTRL+W
Exiting the software to exit software, select File→Close Form or ALT+F4.

Edit Menu

REPORT INFORMATION You can edit referring doctor, examiner, interpreting doctor, and practice information and logo under this menu item. Logos and signatures must be of BMP format to be loaded into the system. There are some samples pre-filled, and you can add, edit, or delete these as you deem necessary.

Add/Delete/Edit Doctor/Examiner/Practice Information/Logos

Under the Edit menu there are multiple options for editing the configuration of your software. You may add, delete, and edit practice information, examiner information, doctor information, output formats, impressions, and ICD codes. You can also upload practice logos to customize reports.



Video Tutorials



[Importing & Converting a Logo](#)

[Reporting Output Formats](#)